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Eluents for Degassers

Even small quantities of other substances, such as additives, modifiers, or salts can influence the durability of the materials.

Note: The list of selected solvents was compiled based on research in the pertinent literature and is only a recommendation by the manufacturer. If there is any doubt, contact the technical support.

Not suitable eluents

- concentrated acids
- concentrated bases
- eluents containing particles
- halogenated hydrocarbons, e.g. Freon®
- Heptane
- Hexane
- Hexafluorisopropanol
- Hydrochloric acid
- Hydro fluoro solvents
- perfluorinated eluents, e.g. Fluorinert® FC-75, FC-40
- perfluorinated polyether, e.g. Fomblin®
- Sodium azide
- Toluene

Less suitable eluents

- Methylene chloride
- diluted phosphoric acid
- slightly volatile eluents

Suitable eluents

- Acetic acid, glacial
- Acetone
- Acetonitrile
- Buffers
- Chlorobenzene
- Chloroform
- Dichloromethane
- Dimethyl sulfoxide (DMSO)
- Ethanol
- Ether
- Ethyl acetate
- Isopropanol
- Methanol, anhydrous
- Methyl ethyl ketone
- Methyl isobutyl ketone
- Methyl tert-butyl ether
- N,N-Dimethylacetamide
- N,N-Dimethylformamide
- n-Propyl alcohol
- NMP (N-Methyl-2-pyrrolidinone)
- Salt
- Tetrahydrofuran (THF)
- Water