K-7400

► **Semimicro Osmometer
User Manual**

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CONTENTS

Using this Manual	6
Conventions in this Manual	6
SOP's in this manual	6
Setting up the Osmometer K-7400	7
Unpacking the Osmometer K-7400	7
Delivery Content	7
Power supply	7
Placing the Osmometer K-7400	7
Description of the Osmometer K-7400	8
Front View	8
Rear View	8
Connection of the Measuring Head	8
Measuring Head Storage	9
Measurement Principle and Process	10
Using the Osmometer K-7400	11
Switching on	11
Internal Software Structure	12
Preparing the Osmometer K-7400	13
Calibration	13
Pre- calibration	13
Measuring calibration	14
Measuring Samples	15
Important Hints for Using	16
Possible Error Messages	17
Self-Freezing of the Sample	17
Non-Freezing of the Sample	18
Using the EuroOsmo [®] 7400 Software	18
Installation	18
Hardware Installation	18
Description of the main screen	19
Maintenance by the User	25
Function Checks of Cooling Capacity	25
Spare Parts and Accessories	26
Technical Data	26
Warranty statement	49
Declaration of Conformity	49
KNAUER Software License Agreement	50

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Using this Manual

Conventions in this Manual


Special Warnings are indicated by the marginal warning sign and printed in bold letters.


The marginal lamp symbol indicates **helpful advice**.

Important Hints are marked by the marginal hand symbol.

Arrows in an outlined form like this example used in block diagrams indicate an automatic program run and change to the next line without the necessity of manual interventions.

Arrows like this, used in block diagrams, indicate that the user is asked to press the corresponding arrow keys

The triangles symbolize the use of corresponding arrow keys.

SOP's in this manual

The **Standard Operating Procedures (SOP)** provided with this manual offer a convenient way of structuring complex tasks in the operation of your Semi Micro Osmometer K7400. They include step-by-step instructions leading the user through all routine tasks during operation. They can be used for documentation purposes and be copied, applied, signed, and filed in order to document the performance of the instrument.

Please operate the instrument and all accessories according to the instructions and SOP's in this manual. This ensures proper results and longevity of your equipment.

Table 1 **SOP's in this manual**

SOP 1	Preparing a Measurement	13
SOP 2	Pre- calibration after measuring head exchange	13
SOP 3	Performing a Calibration	14
SOP 4	Performing a Sample Measurement	15

Setting up the Osmometer K-7400

Unpacking the Osmometer K-7400

Unpack the instrument and check the Osmometer K-7400 and accessories for any damage due to shipping. Make any claims for shipping damage to the transportation company responsible for shipping.

Please check that the semi-micro Osmometer K-7400 is complete, using the Delivery Content List as a guide at page 7. Should any article be missing, please contact our service department.

A complete list of available spare parts and consumables you can find in the section „Spare Parts and Accessories“ on page 26.

Delivery Content

The standard delivery contents of:

Semi-micro Osmometer K-7400

Printer (optionally)

User manual

Power supply cable

RS232 cable

Measuring head with thermistor and stir assembly
(version for either glass or plastic vials)

Pack of 100 measuring vessels

Bottle cleaning solution

Pack 10 ampoules calibration solution 400 mOsm/kg

Pack of cotton swabs

Power supply

Make certain that the correct voltage has been set on the rear panel of the instrument, the power supply is grounded and a corresponding 3-pole power cable is used.

Connect the Osmometer K-7400 to the power supply and switch on the instrument. After switching ON the display presents for a short time information's concerning the instrument version and number.

Placing the Osmometer K-7400

During a measurement no greater temperature changes are allowed.

The air convection holes at both sides of the instrument must remain uncovered.

Otherwise the instrument temperature may increase to values, which disclose a proper function of the Osmometer K-7400. This leads to an automatically shutting down of the display.

To avoid serious damages the reason for the overheating has to be omitted.

Description of the Osmometer K-7400

Front View


Fig. 1 Front view of the Osmometer K-7400

Buttons Function

START/STOP

1. Starts a measurement or calibration run.
2. It interrupts any running measurement!
3. If no run can be started or stopped it causes a jump into the main menu.

VIBRATOR

Activates the stirrer for one second. During a measurement or calibration run the button is not in function.

ARROW KEYS

Cursor navigation buttons for parameter or value selection

Rear View


Fig. 2 Rear view of the Osmometer K-7400

Connection of the Measuring Head

Parts of the measuring head are an inbuilt motor with a stirring wire (vibrator), the thermistor and a cone adapting the measurement vial. The measuring head cable is to be plugged into the connector.

The measuring vial, filled with exact 0.15 ml solution (calibration or measuring) is to fix to the adapter. Two versions of measuring heads are available, for the use of glass vials and plastic vials respectively. The

measuring and calibration procedures are independent on the type of the measuring head.


Fig. 3 Measuring head scheme of the Osmometer K-7400


Fig. 4 Measuring heads of the Osmometer K-7400 for glass (A) and plastic vials (B)

The obvious difference of glass and plastic vials consists in the fixing manner. The plastic vials are simply to push on the cone while the cylindrical glass vials are to push into the adapter sleeve. This gives also the necessity for a different supporter length.

It is not possible to transform a head for glass vials into one for plastic vials or vice versa simply by using or removing the adapter sleeve. Due to the different thermal conduction also internal deviations exist.

Any change of the measuring head type causes the need for a pre-calibration to adapt the instrument to the measuring head (see SOP 2 on page 13).

The measuring head is placed on the housing plate with the measuring vial rising into the cooling compartment. The instrument is now ready for the measurement.

Measuring Head Storage

Generally the measuring head shall be stored with attached measuring vial.

For temporarily put down the measuring head without measuring vial while working with the osmometer use the four feet on top of the instrument.


Fig. 5 Putting down the measuring head without vial

Measurement Principle and Process

The osmotic pressure in a solution is proportional the freezing point depression of that solution. An aqueous solution, with an osmotic pressure corresponding to an ideal 1-molal solution, freezes at $-1.858\text{ }^{\circ}\text{C}$. An aqueous solution with this freezing point has a concentration of 1 Osm/kg.

To perform a measurement first the sample solution is cooled down without stirring. Even pure water can be super-cooled down to about -5 through $-8\text{ }^{\circ}\text{C}$ without freezing. The freezing is initiated by an automatically start of the vibrator at a set temperature. The temperature reaches the freezing point, $0\text{ }^{\circ}\text{C}$ in case of pure water.


Fig. 6 Temperature-time-curve for pure water

The freezing point of a solution is below $0\text{ }^{\circ}\text{C}$. The freezing point depression ΔT is a measure for the osmolality of a solution, which can be read directly from the display given in mOsm/kg.


Fig. 7 Temperature-time-curve for a sample solution

During a measurement the actual temperature is displayed. The instrument detects the maximum temperature reached, stores it, and if crystallization occurred satisfactory, displays the result.

After finishing the measurement the sample vessel is warmed up to 20°C. As this is the stand-by temperature of the cooling cavity, the measuring head can remain there before and after measurement as long as desired.

Using the Osmometer K-7400

Switching on

After switching ON the display presents for a short time information concerning the instrument and the firmware version:

KNAUER	V01.XX
Freezing Point Osmometer	

Now the osmometer is ready for use. In the display the main menu is shown with the result of the last measurement. The cursor is placed on the rhombus field ♦.

>>> 1234 mOsmol <<<
♦ <next sample>

Depending on the actual settings the measurement value can also be shown as a percentage or temperature value.

>>> 12.34 % <<<
♦ <next sample>

>>> -1.234 °C <<<
♦ <next sample>

Using the osmometer the first time neither measuring or calibration values are present. Accordingly the last value will not be displayed.

```

>>> _____ mOsmol <<<
◆ <next sample>
  
```

Internal Software Structure

Starting with the cursor on the rhombus field further menus become accessible in an endless loop by pushing the vertical arrow keys ▲ ▼. Pushing the horizontal arrow keys ► ◀ moves the cursor from field to field. Only those fields are accessible, which are allowed for changing the setting. Other fields like the actual measured temperature will be overleaped.

On any changeable field the possible settings are scrolled using the ▲ ▼ keys. The displayed values will be accepted and stored when leaving the field with one of the ► ◀ keys.


Fig. 8 Example menu sequence of the Osmometer K-7400

Screen **2 TEMP** displays the actual measured temperatures of the sample and the cooling compartment. The values can be changed from [°C] to [°F] with the ▲ ▼ keys.

Screen **3 DATE** enables the setting of date and time. The setting remains active when switching of the instrument.

Screen **4 CALIBRATE** displays the calibration points which are actual in use. 3 calibration points are possible (second row). The first field is fixed on „0000“. At least two calibration points are necessary for any measurement. From this screen the calibration runs will be started, refer **SOP 3** on page 14.

Calibration and measurement must be performed at the same settings. After changing the settings, a new calibration is required.

In screen **5 COOL** the initiation point „freeze“ can be set within the range from -9.0 to -4.0 °C and the minimum temperature „cooling limit“ within the range from -18.0 to -12.0 °C, both in steps of 0.5°.

Screen **6 CALIBRATE-T** enables the pre-calibration to adapt the osmometer to the selected measuring head, refer **SOP 2** on page 13.

Preparing the Osmometer K-7400

The instrument should be switched ON at least 5 minutes before use.

The vibrator key may be pressed to check the stirrer. It will run for a second.

Optionally a printer can be connected to the RS232 socket. The result of each measurement will be printed automatically in the format:

Measured 0400 mOsmol / kg
23.04.15 08:24:35

After a calibration run you get the outprint:

Calibrated 1234 mOsmol / kg
23.04.15 07:36:17

In case of interrupted runs (caused by any error) the screen displays the corresponding error message and you get the printout

(NO) FREEZE
23.04.15 09:06:04

without specifying the individual error.

SOP 1 Preparing a Measurement

1. Place 0.15 ml sample or calibration solution into a clean, dry measurement vial.
2. Put the vial all the way into the adapter. The meniscus of the liquid must be horizontal.
3. Place the measuring head on the instrument in such a way that the vial extends into the cooling cavity.
4. Thus the instrument is ready for calibrating or measuring.

Calibration

Pre- calibration

If you intend to use both types of measuring heads alternatively with the osmometer, it is necessary after each exchange to adapt the instrument to the new measuring head by a pre-calibration procedure as described in SOP 2 on page 13.

SOP 2 Pre- calibration after measuring head exchange

1. Prepare the instrument according SOP 1 for a measurement of deionised water.
2. Select the CALIBRATE-T screen.
3. Press the start key. While the pre-calibration is running the screen displays the sample temperature and the elapsed time.

<p>0000 mOsmol > +5.43°C 01:23</p>

At any time the run can be interrupted pressing the STOP key.

4. After finishing the run on the screen is displayed „successful“. The previous measuring calibration will be aborted.

Any performed pre-calibration deletes the existing measuring calibration, because it is not valid for the new measuring head!

5. If you get „not successful“ the sample is not frozen. Please check that
 - the sample is a 0 mOsm/kg one
 - a valid measuring head is proper connected to the instrument and repeat the pre-calibration.

Measuring calibration

Normally a two point calibration is sufficient.

Calibration must be performed very carefully. A confusion of the calibration solutions will affect all later measurements and give incorrect results.

Up to three calibration points are possible, two are necessary. On any of the second row fields the desired calibration concentration can be selected. After this the calibration run is started by pressing the START/STOP button.

The first field is fixed on „0000“. On the other fields the choices are ----, %%%%, 50, 100, 150, up to 1900, 1950, 2000 in steps of 50 mOsmol/kg..

Selecting „----“ the corresponding calibration point will be deactivated.

The selection of „%%%%“ deactivates all other calibration points except „0000“. Simultaneously the displayed measurement results will be shown in % of this calibration point instead of „mOsmol/kg“.

At least two calibration points must be active, the „0000“ and one of the others.

The calibration point remains active until a new valid value for the selected field is created. If a calibration run is not successfully finished the last calibration value becomes reactivated.

SOP 3 Performing a Calibration

1. Prepare the instrument according **SOP 1** for a measurement of deionised water.
2. Set the first field in the second row of the CALIBRATE screen to „0000“.
3. Press the start key. While the calibration is running the screen displays the selected calibration point, the sample temperature and the elapsed time.

0000 mOsmol > +5.43°C 01:23

At any time the run can be interrupted pressing the STOP key.

4. After finishing the run on the screen is displayed

0000 mOsmol > -0.00°C< o.k.? <start>

To accept this value presses the START key again. An earlier measured calibration will be overwritten. Pushing any other key the

new value will be abandoned and the CALIBRATE menu will be shown again.

5. Wait for the warming up to the stand-by temperature.
6. Prepare the instrument according **SOP 1** for a measurement of a 400 mOsm/kg calibration solution (KNAUER Order No. Y1241).
7. Activate the second calibration field and select "0400" and repeat steps 3 to 5.

If for special applications as some times in the pharmaceutical industry a 3 point calibration is required proceed according **SOP 3** and repeat the procedure for calibration field 3.

A new calibration is also necessary after

1. exchanging the measurement head
2. changing the settings of FREEZING potentiometer

Measuring Samples

SOP 4 Performing a Sample Measurement

1. Prepare the instrument according SOP 1 for a measurement of a sample solution.
2. Press the START key to start the measurement run. While running the screen displays the sample temperature and the elapsed time.

RUN >	-1.23°C<
01:23	

At any time the run can be interrupted pressing the STOP key.

3. The instrument detects the freezing point depression and if the crystallisation occurred properly displays the corresponding osmolality. On the main screen.

>>> 1234 mOsmol <<<
◆ <next sample>

4. The printer output is given in the format

Measured 1234 mOsmol/kg
19.03.15 16:05:35

If no crystallisation (temperature increase) takes place after the automatically start of the vibrator the measurement will be stopped with the error message

> ERROR <
NO FREEZE

If the crystallisation (temperature increase) takes place before the start of the vibrator the measurement also will be stopped with an error message

> ERROR <
FREEZE

Do not remove the measurement vessel from the adapter before the sample is molten. Otherwise the danger of damaging the measuring head is given.

The sample melting can be accelerated by warming up with the fingers.

Important Hints for Using

1. To ensure a fixed sample volume, always introduce the solutions with a clean and dry pipette tip into the measuring vessel.
2. Handle the thermistor always very carefully. All abrasive materials should be kept away from it. Measuring vessel and thermistor should be cleaned occasionally with KNAUER cleaning solution supplied with the instrument.
3. (Dilution: 1:10 solution in distilled water of 40 ± 5 °C)
4. A dirty thermistor does not give accurate measurement results!
5. Cleaning and drying the measuring vessels with alcohol or acetone is not recommended because remaining solvent vapours could cause false results!
6. Solutions containing proteins, such as sera, can only be measured once. Freezing causes denaturation of the protein so that repetition of the measurement with the same sample would result in increased values. Thus, in the case of sera, only one measurement is possible for one sample.
7. Please note the difference between deal and real osmolality. Only real osmolality can be measured. It is not possible to prepare other standard solutions by dilution of a calibration solution since the activity coefficient of the solution changes with dilution.
8. It is sometimes difficult to cool distilled water which is kept in polyethylene bottles for a longer period. We recommend in this case storing the water which is used for the zero point calibration in glass bottles.
9. If the same salt solution is measured several times, it must be mixed after thawing (stir briefly). During the thawing process ice floats to the top of the solution. As ice does not contain salt, the top layer of the solution is diluted while thawing.
10. Calibration solutions can become more concentrated if the bottle is opened more than 1 hour. We recommend using KNAUER calibration solutions in glass ampoules (Order No. Y240, Y1241, Y1248, Y1250)
11. Depending on the room conditions it may be possible that after a longer period in the standby mode a little bit of water is condensed in the cooling compartment. Prior to the next measuring it should be removed with a dry cloth.

Possible Error Messages

Any error message will be deleted simply by pressing any key. The main screen will be displayed.

In any case the first row shows > ERROR < which is specified in the second row.

Table 2

Message	Cause	Removal
FREEZE	The crystallisation took place or another fault occurred before the start of the vibrator	See Table 3 On page 16
NO FREEZE	No crystallisation took place after the automatically start of the vibrator	See Table 4 on page 18
OUT OF RANGE	The measured value is outside the acceptable range	Usually malfunction of the thermistor
NOT CALIBRATED	No calibration points are available	See SOP 3 on page 14
COOLER TIMEOUT	Timeout due to too slow cooling.: >3 minutes	
MEASUREMENT TIMEOUT	Timeout due to too long measurement.: >5 minutes from START	

Self-Freezing of the Sample

If the sample does not super-cool but freezes continuously without any external agent there are three possible causes:

Table 3

Error reason	Solution
Too less sample solution is in the measuring vial. A very cold zone, which causes self-freezing of the sample, develops above the solution on the measuring vial material.	Use exact the required sample volume of 0.15 ml
Scratches or contaminations at the thermistor or measuring vial nucleate inducing premature freezing.	Exchange of the measuring vial and cleaning the thermistor with cleaning solution.
Solutions saturated with air or gases are sometimes difficult to super-cool.	Degas the sample solution by ultra sound

If these three remedies do not produce the desired result, it is indicated that the difficulties are caused by the nature of the sample to be measured and the set cooling temperature is too low. In this case the freezing temperatures must be set higher. Setting the cooling temperature see screen 5 COOL in Fig. 8 on page 12.

Calibration and measurement must be performed at the same setting. After changing the setting, a new calibration is required.

Non-Freezing of the Sample

If the sample does not freeze, one of the following causes may be responsible:

Table 4

Error reason	Solution
The cooling rate is not sufficient	Ensure that the air inlets and outlets at the instrument are free.
Too much sample solution is in the measuring vessel leads to a no cooled zone in the upper part of the solution. The stirring results only in mixing of warm and cold sample layers.	Use exact the required sample volume of 0.15 ml

If these remedies do not produce the desired result, it is indicated that the difficulties are caused by the nature of the sample to be measured and the set cooling temperature is too high. In this case the freezing temperatures must be set lower. Setting the cooling temperature see screen 5 COOL in Fig. 8 on page 12.

Calibration and measurement must be performed at the same setting. After changing the setting, a new calibration is required.

Using the EuroOsmo® 7400 Software

The whole features of the osmometer K-7400 become available, running it with the software EuroOsmo® 7400 for Windows, version 1.5 or higher.

The License Agreement should be read before using the EuroOsmo® 7000 software! Using the software, the user indicates that he has read the License Agreement and agrees with the regulations.

This chapter will explain the basic knowledge needed for working with K-7400 and the EuroOsmo® 7400 software.

Installation

The hardware requirements (processor, working memory) are very low for the installation of EuroOsmo® 7400. The program runs under Windows7 and 8 but also with older operating systems, back to Windows XP.

1. Insert CD into CD drive of the PC
2. Wait up to 30 seconds. The set-up starts automatically
3. If your CD ROM drive does not have AutoPlay:
 - 3.1. Double-click on „My Computer“
 - 3.2. Right-mouse click on your CD ROM
 - 3.3. Select AutoPlay from the right-click menu
 - 3.4. If there is no AutoPlay menu, select Explore from the right-click menu, go to „EuroOsmo 7000/Disk1“ and double-click on „Setup.exe“
4. Follow the installation guide.

Hardware Installation

The Osmometer K-7400 has to be connected with a defined port of the computer (com1 for instance) prior to starting the program.

Description of the main screen

After starting the program usually by a double-click on the corresponding icon you will get a screen similar to Fig. 9.


Fig. 9 Opening screen of the EuroOsmo® 7400 Software

If the osmometer is not connected or switched off you will get the corresponding message.


Fig. 10 Error message for a not connected osmometer

Table 5 Description Screen Elements

Menu line	File Edit View Help
File	Menu for opening, saving, and printing of data files or previews (measurement or calibration)
Edit	Menu for editing file changes and search performing.
View	Menu to show or to hide status bar and toolbar
Help	The Version of the EuroOsmo® is shown
Toolbar	
	Most of the icons are usual windows features and do not need further explanation
	Opens in any case a printer setup window prior printing
	Displays a preview of the actual report
	Starts the search function of EuroOsmo 7400
	Paste the data of a selected (first) field into following lines of the chosen column
	After each bar code scan the cursor switches to the next sample line

Run Settings

Sample Number:

Select the number of the sample to be measured via the pull down menu.

Start button for any measurement or calibration run. While starting the run it changes to a stop button to enable a manually aborting of the run.

If the sample has already been measured the following warning will be displayed when the button is used.



Fig. 11 Warning bevor overwriting existing results

Sample List

This list is used to name the samples and their codes. The measurement results will be filled in automatically. Using the scroll bar at the bottom a **Comments** column becomes accessible for detailed information. Using the „insert row“ and „delete row“ functions of the EDIT menu the list can be extended ore shortened.

The samples will be measured in the order of the list and the results are filled in and stored in it.

The screenshot shows a window titled "Sample List" with a menu bar containing "Sample List", "Setup", and "Report & Graph". The "Sample List" menu item is circled in red. Below the menu bar is a table with the following data:

No.	Sample Name	Sample Code	Result	Units
1	Sample 1	Code 1		
2	Sample 2	Code 2		
3	Sample 2	Code 2		
4	Sample 2	Code 2		
5	Sample 3	Code 3		
6	Sample 4	Code 4		
7	Sample 5	Code 5		
8	Sample 8	Code 8		
9	Sample 9	Code 9		
10	Sample 10	Code 10		

Fig. 12 New sample list

Double clicking on the No. of an already measured sample displays the corresponding temperature sequence. On top of the graph the path of the already stored sample list is shown.


Fig. 13 Measuring window with the temperature sequence of a run


Fig. 14 Setup area of the main window

Setup Area for entering general settings for measurement sequences and corresponding reports.

Operator Enter the operator's name

Communication Port Select the COM port to which your osmometer is connected.

Working Parameters

Monitoring view The entered value defines the scale of the monitor by its final value in minutes. The temperature pattern is shown continuously for the defined period.

Measurement view The entered value defines the scale of the measurement window by its final value in minutes.

Freezing Level The initiation point „freeze“ can be set within the range from -13.0 to -4.0 °C. Any decimal value will be rounded in steps of 0.5°.

Cooling Limit The minimum temperature can be set within the range from -18.0 to -12.0 °C. Any decimal value will be rounded in steps of 0.5°.

Units Select the units valid for the measurement sequence. The choice is mOsmol, °C, °F, and %. If °F is selected the values of the freezing level and cooling limit are automatically changed to °F.

By clicking on the Apply button the entered values are sent to the osmometer. If any value is out of the allowed range a corresponding message will

be displayed. The setting of the osmometer remains unchanged.

Fig. 15 Error message in case of not valid parameters

If a project file *.osm was loaded the stored settings will be displayed but not sent to the osmometer. To achieve this press the **Apply** button.


Fig. 16 Report design area

Report & Graph

Report Settings Report design area

Logo Bitmap Clicking on the  button an explorer window for the selection of the file with the logo to be used in the report.

Change Caption Clicking on the  button a window for the report caption design will be opened.


Fig. 17 Setup window for the report caption

Enter the wording of the report caption and click on the **Font...** button to format the caption.



Fig. 18 Setup window for the report caption format

Change SubTitle Clicking on the  button a window for the report sub title design will be opened, analogue to Fig. 17.

Graph Settings Area for settings of the graph design.

A click on one of the four buttons will open the colour selection window (Fig. 19) where the colours of the curve, axes, grid lines, and the background can be defined separately.

Fig. 19 Color selection window

Grid On This option field is to show or to hide the grid lines of the graphs.

The click on the  icon opens a printer setup window prior printing the report. A preview of the actual report will be displayed when clicking on the  icon.

Search function

The file super ordinate search function enables the locating of individual samples using the operator, sample name, sample code, date or a value range of results.


Fig. 20 Search function window

Clicking on the  button an explorer window for the selection of the path for the search will be opened. Additionally you have the option to look also in subfolders or not. Search subjects are the operator, the sample name, and/or the sample code. As wildcard characters can be used for single characters and * for not defined areas before and/or following to the defined search characters.

If necessary the search results can be reduced by advanced options. These are the measurement values and the date of the measurement.


Fig. 21 Search function window, advanced options

Maintenance by the User

Function Checks of Cooling Capacity

The cooling capacity can be checked by measuring the temperature in the cooling cavity (-13 to -16°C during cooling down and ~ 20°C in stand-by position).

To perform this test, the instrument must be calibrated. Connect the measuring head cable and place it on top of the instrument (Fig. 5). Press Start/Stop to start the measurement and wait 5 minutes to get the error message <MEASUREMENT TIMEOUT>. After this press quickly twice the ∇ key to read out the *cooler* value on *cool* menu (Fig. 8.2).

The value must be in the range of -13 to -16°C. If the value is not reached, please check that the *cooling limit* is set to -16°C (Fig. 8.5). If you use for this test a different *cooling limit* than -16°C you get different results. Retry the test and press quickly twice the ∇ key after 5 minutes without displaying the error message. If you change the *cooling limit*, you must perform a new calibration after this test.

Additionally a calibrated reference thermometer can be used for measuring the temperature in the cooling block of the osmometer. Please verify, that the reference thermometer has contact to the cooling unit and displays quickly temperature changes. Fill the copper cooling cavity close to the top with pure water and install the thermometer. Perform the test as described above and observe the reference thermometer. After the test, it is necessary to dry completely the copper cooling room with a tissue.

The cooling capacity is usually sufficient for a sample volume of 0.15 ml if the freezing turbulence in the solution is initialized by a start of the vibrator within 150 seconds.

If the temperature of -13°C cannot be obtained, the Osmometer K-7400 must be returned for repair.

The thermistor should always be placed in the sample centre 0.5 to 1 mm above the bottom of the measurement vial.

Be sure not to bend the vibration wire.

Fill 150 μ l of a sample or calibration solution into the measurement vial. Carefully insert the measuring vial all the way into the adapter. The thermistor should now be approximately 0.5 to 1 mm above the bottom of the measurement vessel.


Fig. 22 Thermistor Adjustment

Spare Parts and Accessories

	KNAUER Order No.
Power cable, 230V	M1479
RS232 connection cable	M1702
Manual	V7024
Plain printer	A3711
Measuring head for measuring glass vials	A0840
Measuring head for measuring plastic vials	A0840-1
Adapter for disposable measuring glass vials	A0915
Pack of 100 measuring glass vials	A0914
Pack of 1000 measuring glass vials	A0913
Pack of 500 measuring plastic vials	A0272
Pack of 1000 measuring plastic vials	A0720
Pack of 12 vials per 5 ml calibration solution 300 mOsm/kg	Y1240
Pack of 10 vials per 2 ml calibration solution 400 mOsm/kg	Y1241
Pack of 12 vials per 5 ml calibration solution 850 mOsm/kg	Y1250
Pack of 12 vials per 5 ml calibration solution 2000 mOsm/kg	Y1248
Pack of cotton swabs	A0469
Paper roll for plain printer	V7013

Technical Data

Sample volume:	0.15 ml
Solvent:	water
Measurement time:	approx. 2 min
Ambient temperature:	40°C maximum
Display:	LCD, 2 lines 24-characters each
Accuracy:	≤ 1% (for standardized measuring vessels)
Precision:	RSD < 1% (for standardized measuring vessels)
Measuring range:	0 – 1000 (2000) mOsm/kg
Calibration:	2 point calibration (0 and one free selectable osmolality)
optional:	3 point calibration (0 and 2 free selectable osmolalities)
Measuring procedure:	Completely automatic after sample insertion
Power supply:	90-260 V, 47 – 63 Hz, 70 VA
Fuses:	T 1 A, DIN 41571
Dimensions:	160 x 182 x 340 mm (WHD)
Weight:	5.0 kg

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Warranty statement

The factory warranty for the device is stipulated by contract. This warranty is valid for one year from the date of shipment. Warranty is void by improper operation, application and/or maintenance of the instrument. Any expendable components such as glass parts and thermistors are excluded from this warranty. This warranty is validated when, within fourteen days of receipt, the accompanying guarantee card is completed and returned to the manufacturer:

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If we find a defect covered by the guarantee, repair or replacement, at our discretion, will be carried out free of charge. Packing and transport costs are borne by the purchaser.

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Declaration of Conformity

The Declaration of Conformity is part of the delivery and accompanies the product as a separate document.

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KNAUER Software License Agreement

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(hereinafter Licensor or KNAUER)

and the software user (hereinafter Licensee)

for the software package

EuroOsmo[®]7400 V1.5 and handbooks.

In return for acquiring a license to use the software ("EuroOsmo[®]7400 V1.5") and related documentation, Licensee agrees to the following terms and conditions:

1. License.

This agreement grants the Licensee a license to

- (a) use the EuroOsmo[®]7400 V1.5 on a single computer system or, in the case of a multi-user or networked system which permits access to the EuroOsmo[®]7400 V1.5 by more than one user at the same time, at one (1) single working location, and
- (b) make one copy of the EuroOsmo[®]2000 V1.2 in machine readable form solely for back-up purposes provided you reproduce KNAUER's copyright proprietary legends and mark the back-up copy with the express wording "back-up copy Knauer GmbH, Hegauer Weg 38, D-14163 Berlin".

2. Restrictions.

Licensee may not distribute copies of the EuroOsmo[®]7400 V1.5 to others or electronically transfer the EuroOsmo[®]7400 V1.5 from one computer to another over a network. Licensee may not use the EuroOsmo[®]7400 V1.5 from multiple locations of a multi-user or networked system at any time. The EuroOsmo[®]7400 V1.5 contains trade secrets and in order to protect them the Licensee may not decompile, reverse engineer, disassemble, or otherwise reduce the Software to a human-perceivable form: LICENSEE MAY NOT MODIFY, ADAPT, TRANSLATE, RENT, LEASE, LOAN, RESELL FOR PROFIT, DISTRIBUTE, NETWORK OR ANY DERIVATIVE WORKS BASED UPON THE SOFTWARE OR ANY PART THEREOF. The Licensee must not remove from the program firm names, trademarks, logos, copyright notices and other markings or devices, which serve the protection of the program.

3. Ownership of EuroOsmo[®]7400 V1.5.

Licensee owns the media upon which the EuroOsmo[®]7400 V1.5 is recorded or fixed, but KNAUER, Licensor, retains title and ownership of the EuroOsmo[®]7400 V1.5 recorded on the original media and all subsequent copies of the EuroOsmo[®]7400 V1.5, regardless of the form or media in which or on which the original and other copies may exist. This license is not a sale of the EuroOsmo[®]7400 V1.5 or any copy.

4. Confidentiality.

Licensee agrees to maintain the EuroOsmo[®]7400 V1.5 in confidence and to not disclose the EuroOsmo[®]7400 V1.5 to any third party without the express written consent of KNAUER. Licensee further agrees to take all reasonable precautions to preclude access of unauthorized persons to the EuroOsmo[®]7400 V1.5. The Licensee must not pass copies of the program or the handbook to third persons.

5. Term.

This license is effective until January 1, 2030, unless terminated earlier. Licensee may terminate the license at any time by destroying the

Software (including the related documentation) together with all copies or modifications in any form. KNAUER will have the right to terminate the Licensee's license immediately if the Licensee fails to comply with all copies or modifications in any form. The termination does not raise Licensee's claim for refunding of the license fee.

6. Limited Warranty.

KNAUER warrants only that the media upon which the EuroOsmo[®]7400 V1.5 is furnished will be free from defects in material or workmanship under normal use at the time of delivery. KNAUER DOES NOT AND CANNOT WARRANT THE PERFORMANCE OR RESULTS LICENSEE MAY OBTAIN BY USING THE SOFTWARE OR DOCUMENTATION. THE FOREGOING STATES THE SOLE AND EXCLUSIVE REMEDIES KNAUER WILL PROVIDE FOR BREACH OF WARRANTY. EXCEPT FOR THE FOREGOING LIMITED WARRANTY, KNAUER MAKES NO WARRANTIES; EXPRESSED OR IMPLIED, AS TO NON-INFRINGEMENT OF THIRD PARTY RIGHTS, MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. Some countries do not allow the exclusion of implied warranties or limitations on how long an implied warranty may last, so the above limitations may not apply to the Licensee. This warranty gives the Licensee specific legal rights and the Licensee may also have other rights which vary from country to country.

7. Limitation of Liability.

IN NO EVENT WILL KNAUER BE LIABLE TO LICENSEE FOR ANY SPECIAL DAMAGES, INCLUDING ANY LOST PROFITS, LOST SAVINGS, OR OTHER INCIDENTAL OR CONSEQUENTIAL DAMAGES - EVEN IF KNAUER HAS BEEN ADVISED OF THE POSSIBILITY OF SUCH DAMAGES, OR FOR ANY CLAIM BY ANY OTHER PARTY. Some countries do not allow the exclusion of limitation of special, incidental, or consequential damages, so the above limitation, or exclusion may not apply to the Licensee. Any liability is – in any case – limited to the license fee.

8. Limitation of Remedies.

KNAUER's entire liability and Licensee's exclusive remedy shall be:

- (a) the replacement of any media not meeting KNAUER's limited warranty which is returned to KNAUER, or
- (b) if KNAUER or its distributor is unable to deliver replacement media which is free of defects in materials or workmanship, the Licensee may terminate this Agreement by returning the EuroOsmo[®]7400 V1.5 and the Licensee's money will be refunded.

9. This license is not assignable.

The Licensee does not have the right to assign rights or grant sublicenses.

10. License Fee / Royalty.

For the license to use the contractual EuroOsmo[®]7400 V1.5 in the above defined scope and within the restrictions the Licensee pays to the Licensor a fee. The amount becomes due and payable upon delivery of the program discs. The license to use of the program commences with payment of the royalty.

11. If a provision of this Agreement should be invalid or become invalid, then the legal effect of the other provisions shall remain unchanged. Instead of an invalid provision, a valid provision shall be deemed to be agreed which comes closest to what the parties intended economically. This document is the complete Agreement between KNAUER and Licensee; it supersedes all and any prior agreement or offer. All and any modifications of this agreement have to be made in writing. For all controversies out of this Agreement, the District Court Berlin shall have exclusive jurisdiction. The applicable law shall be that of the Federal

Republic of Germany (German Civil Code) excluding the UN-Code on the international sale of goods.