HPLC Pump 1800 Manual / Handbuch

V 7115 10/2008

Smartline

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CONTENTS

General Description	3
Using this Manual	3
Conventions in this Manual	3
SOP's in this manual	4
Setup of the Preparative Pump 1800	5
Unpacking	5
Front View of the Preparative Pump 1800	5
Identifying installed Pump Heads	6
Pump Heads and Hydraulic Connections	7
Rear Panel View of the Preparative Pump 1800	8
Connection to Main Power	8
Control of the Preparative Pump 1800	9
Basics for operating the Preparative Pump 1800	9
Power On and Self-test	9
Appearance of the Cursor	9
Menu Structure of the Preparative Pump 1800	9
Using the Arrow keys in the SETUP Menu	10
Using the SETUP Menu	12
Using the VIEW Menu	13
Using the GLP Menu	14
Using the LINK Menu	15
Operating the Preparative Pump 1800	15
Working Example for setting parameters	15
Using the Low Pressure Gradient Mode	16
Using the High Pressure Gradient Mode	18
Direct Control of the Pump	19
Maximum Pressure	19
Events	20
Automatic Piston Backflushing	21
Programming the Preparative Pump 1800	22
Features of Programs	22
Example for creating a program	25
Using Program Links	26
Connecting the Preparative Pump 1800 to Other Instruments	27
Hydraulic Connections	27
100, 250 and 500 ml Head	27
1000 ml Head	27
Electrical Connections	28
Using the EVENT and ANALOG/ERROR/START Sockets	28
Configuring Plug Strips	28
Analog control of the pump	30
Using the RS232 Serial Interfaces in the KNAUER-Net	30
Easy Maintenance	32
Changing and Cleaning the Pump Head	32
List of Display Messages	38
List of Error Codes	38
Packing List	39
Spare Parts and Accessories	40
Technical Data	41
Warranty statement	82
Declaration of conformity	83
INDEX (english)	84

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General Description

The Preparative Pump 1800 features a double-piston design with active pulsation dampening. It is capable of delivering solvents with very low residual pulsation. It can be configured with three different pump heads to meet a broad range of requirements of flow ranges:

- 100 ml pump head (0.1 – 99.9 ml/min), stainless steel, up to 400 bar
- 250 ml pump head (0.1 – 249.9 ml/min), stainless steel, up to 200 bar
- 500 ml pump head (0.1 – 499.9 ml/min), stainless steel, up to 100 bar
- 1000 ml pump head (0.1 – 999.9 ml/min), stainless steel, up to 50 bar

Pump heads can be interchanged by loosening only four screws. Experienced users can disassemble the heads for simple maintenance tasks like cleaning, changing piston seals, etc.

The easy-to-learn user interface of the Preparative Pump 1800 provides access to the integrated software. It supports quaternary high pressure gradient systems via a master/slave configuration, where up to three additional Preparative Pumps 1800 can be controlled.

Most KNAUER Software Packages, like EuroChrom® 2000, allow for digital remote control of any configuration of KNAUER solvent delivery systems, enabling full and transparent integration of the Preparative Pump 1800 pump into a wide range of chromatographic systems.

A quaternary low pressure gradient former is available.

Using this Manual

Conventions in this Manual


Special Warnings are indicated by the marginal warning sign and printed in bold letters.

The marginal lamp symbol indicates **helpful advice's**.

Important Hints are marked by the marginal hand symbol.

Arrows in an outlined form like this example used in block diagrams, indicate an automatic program run and change to the next line without the necessity of manual interventions.

Arrows like this, used in block diagrams, indicate that the user is asked to press the corresponding arrow keys

The triangles symbolize the use of corresponding arrow keys.

SOP's in this manual

The **Standard Operating Procedures (SOP)** provided with this manual offer a convenient way of structuring complex tasks in the operation of your Preparative Pump 1800. They include step-by-step instructions leading the user through all routine tasks during operation. They can be used for documentation purposes and be copied, applied signed, and filed in order to document the performance of the instrument.

Please operate the instrument and all accessories according to instructions and SOP's in this manual. This ensures proper results and longevity of your equipment.

Table 1 *SOP's in this manual*

SOP 1	Selecting Menus	10
SOP 2	Entering Data	10
SOP 3	Selection of correct pump head:	16
SOP 4	Selection of operation mode	16
SOP 5	Installing the ternary LPG valve block	17
SOP 6	Installing the binary LPG valve block	18
SOP 7	Setting the Maximum Pressure	19
SOP 8	Setting the Flow Rate	19
SOP 9	Setting the Solvent Composition	20
SOP 10	Setting Events	20
SOP 11	Start and Stop of Solvent Delivery	21
SOP 12	Purging the Preparative Pump 1800	21
SOP 13	Preparing your Preparative Pump 1800 for automatic piston backflushing	21
SOP 14	Creating a Program	23
SOP 15	Execution of Programs	23
SOP 16	Modifying and Deleting Programs	24
SOP 17	Using a Wake-Up Program	24
SOP 18	Creating Program Links	26
SOP 19	Executing Program Links	27
SOP 20	Mounting Plug Strips	28
SOP 21	Removing the Pump Head	32
SOP 22	Removing and Checking Piston Rods	32
SOP 23	Disassembling the Pump Head I	33
SOP 24	Disassembling the Pump Head II	34
SOP 25	Assembling the Pump Head and Exchanging Piston Seals I	35
SOP 26	Assembling the Pump Head and Exchanging Piston Seals II	35
SOP 27	Installing the Pump Head	36
SOP 28	Cleaning and Replacing Check Valves	36

Setup of the Preparative Pump 1800

Unpacking

Unpack the instrument and check pump and accessories for any damage due to shipping. Place any claims referring to damage at the transportation company responsible for shipping.

Please check that the Preparative Pump 1800 is complete, use section „Packing List“ on page 39. Should, in spite of our thorough control procedures, something be missing, please contact our service department.

Please remove the transparent protective foil from the pump outlet and the display.

For the complete list of spare parts available please refer to chapter „Spare Parts and Accessories“ on page 40.

Front View of the Preparative Pump 1800


Fig. 1 Front View of the Preparative Pump 1800

The display shows information about the current status of the pump, such as the flow rate, the measured pressure or the current status of the active program. It serves as a control instrument when entering data.

Pumping can be started or stopped by pressing the key START/STOP, see {2.1}, using the preset values. The foil key PURGE, see {2.2}, allows immediate access to the purging function of the pump. Purging has to be stopped manually by pressing PURGE again or by pressing the START/STOP key.

You can use the numeric keys for entering data and control the execution of programs. Detailed information is given in chapter „Control of the Preparative Pump 1800“ on page 9 and „Programming the Preparative Pump 1800“ on page 22.


Fig. 2 Display information and Foil Key Area

Identifying installed Pump Heads

Maximum flow rate, 100, 250, 500, or 1000 ml/min, is indicated on a label on the front of the pump head.

Fig. 3 Pump head label (1000ml head)

The pump heads with inlay design bear in stainless steel version on both inlays only the imprint 100, 250, 500, or 1000. The imprint informs on the total capacity of the pump head. They are not allowed to be added.


Fig. 4 Pump head identification (100ml head, inlay design)

The biocompatible titanium version as additionally signed with Ti.

Pump Heads and Hydraulic Connections


All hydraulic connections are located on the front panel of the Preparative Pump 1800. For connecting eluents to the pump head see section „Connecting the Preparative Pump 1800 to other Instruments“ Fig. 24 and Fig. 25 on page 27.

- | 3 Pump Head | 4 Transducer with Purge Valve |
|---|--|
| 3.1 Connection capillaries | 4.1 Deaeration capillary, (Purging outlet) |
| 3.2 Piston backflushing connectors | 4.2 Deaeration screw for opening purge valve |
| 3.3 Eluent outlets to transducer | 4.3 Eluent outlet to column |
| 3.4 Outlets valve housing | 4.4 Inlets to transducer |
| 3.5 Head set screws | |
| 3.6 Labels indicating the maximum flow rate | |
| 3.7 Backflushing pump inlet | |
| 3.8 Eluent inlets and valve housing | |
| 3.9 Backflushing pump outlet | |

Fig. 5 Pump Head and transducer with Purge valve

Rear Panel View of the Preparative Pump 1800

5 Rear panel elements

- | | |
|------------------------------------|---|
| 5.1 Ventilator | 5.6 Power connector line in |
| 5.2 Remote terminal strip | 5.7 Fuse compartments |
| 5.3 Socket for Ethernet connection | 5.8 Power connector line out |
| 5.4 RS232 Interfaces | 5.9 Connector for optional Low Pressure Gradient unit |
| 5.5 Event terminal strip | 5.10 Serial Number |

Fig. 6 Rear Panel View of the Preparative Pump 1800

Connection to Main Power

The Preparative Pump 1800 is equipped with an universal power supply for input voltage ranges of 90 - 130 Volts AC and 170 - 260 Volts with 40-60 Hz AC. A manual setting of the supply voltage is not required.

Make sure that the main power supply is properly grounded and a correct three-pronged power cable is used. Connecting the instrument to faulty power sockets can cause damages.

Connect the power cable to the socket on the rear of the instrument, see {5.6} in Fig. 6, and switch on the device with the „ON / OFF switch“, Pos. {6} at the front side of the device (Fig. 1).

Further devices with a maximum power consumption up to 500 Watts may be connected to the „Power connector - line out“, Pos. {5.8} in Fig. 6.

Control of the Preparative Pump 1800

Basics for operating the Preparative Pump 1800

Power On and Self-test

On power on the instrument displays first the VERSION of the kernel, followed by the VERSION of the software, e.g. 2.xx. The instrument performs a self test, including firmware validation and tests of RAM and motor unit. After completion the message **SYSTEMTEST OK** indicates that the pump is ready for operation and the ID of the last active program is displayed.

Never run the Preparative Pump 1800 without liquid in the pump head or in the piston backflushing compartment. For piston backflushing, follow the instructions given in chapter „Piston Backflushing“ on page 21. Operating the pump without solvents may lead to damage of the pump seals.

Appearance of the Cursor

You can use the yellow arrow keys on the foil keypad for positioning the blinking cursor on any field of the menu screen.

The cursor appears as underscore character while entering data, see Fig. 7. Use „Numeric Key pad for data input“, Pos. {2.3} in Fig. 2, and refer to SOP 2 „Entering Data“ on page 10.

Fig. 7

Blinking Cursor

Cursor during data input

Menu Structure of the Preparative Pump 1800

The display of the Preparative Pump 1800 uses one main screen with status information of the instrument. From the main screen you can access four additional screens for display of parameters.

By default, the Preparative Pump 1800 displays the parameters of the first program line of the currently active program.

The center part of Fig. 8 „Menu Structure of the Preparative Pump 1800“ shows an example of the main screen.


Fig. 8

Menu structure of the Preparative Pump 1800

SOP 1 Selecting Menus

This SOP applies to the Preparative Pump 1800 Firmware Revision 2.xx. Use this SOP for selecting menus in your Preparative Pump 1800.

You can access the four menus of the Preparative Pump 1800 with the following procedure:

1. Place the cursor in a field in one of the four corners of the main screen, e.g. on **OFF** in the lower left corner of the main menu.
2. Depending on the selected corner, press the corresponding left or right arrow key longer than 2 seconds. In the example above selection of key \leftarrow will activate the VIEW Menu on the display.
3. Select the rhombus \blacklozenge in the lower left corner of the display if you want to access additional pages of the selected menu. You can access all menu pages by scrolling using ∇ or \blacktriangle arrow keys.
4. Select the rhombus \blacklozenge and use \triangleright or \triangleleft to return to the main screen.


Fig. 9 Example of a multi-page menu: VIEW menu

Menus without rhombus are single page menus.

SOP 2 Entering Data

This SOP applies to the Preparative Pump 1800 Firmware Revision 2.xx. Use this SOP for entering parameters and program steps in your Preparative Pump 1800.

1. Place the cursor at the desired position of the main screen.
2. Press a numeric key to start data input. The appearance of the cursor changes to a underline character, see Fig. 7 on page 9.
3. Press numeric keys for entering desired numeric values. Wrong entries can simply be overwritten by continuing the entry.
4. Confirm value by pressing any of the arrow keys. This finishes data entry mode and moves the cursor to the next field according to the arrow key direction.
5. You can modify and delete data according to SOP 16 „Modifying and Deleting Programs“ on page 24.

Using the Arrow keys in the SETUP Menu

Fig. 10 „Structure of the SETUP Menu“ on page 11 gives a complete overview of the Setup Menu.

- You can use \triangleright for accessing menus.
- You can use \blacktriangle and ∇ for scrolling through the list of parameters of any individual menu item.


Fig. 10 Structure of the SETUP menu

Using the SETUP Menu

The SETUP menu of the Preparative Pump 1800 enables access to the following features and functions:

GRADIENT MODE: The Preparative Pump 1800 is capable of isocratic (none), high pressure gradient (**HPG**) and two kinds of low pressure gradient operations (**LPG→ 220** or **LPG→ 1000**).

To select the appropriate gradient mode proceed according to SOP 4 on page 16.

Please observe that both the gradient blocks A4094 and A4099 need the corresponding settings on the pump (LPG→ 220 or LPG→ 1000). Details will be described in the section Using the Low Pressure Gradient Mode on page 16.

VALVES: This window is only for information purposes. The settings can not be changed. The values inform you, which valve block can be controlled. These settings are without any influence while working in isocratic or high pressure gradient mode. If you have changed the LPG→ 220 or LPG→ 1000 setting the corresponding valve values become active with a pump reset (new start).

ternary block A4099				binary block A4094			
VALVES:	Wait	Boost	Hold	VALVES:	Wait	Boost	Hold
◆	10ms	01.5ms	045%	◆	10ms	10.0ms	080%

Fig. 11 Information display for valve control

PUMP HEAD: You can select the type of pump head installed in your pump. Selection of any pump head sets specific flow rate calculating data as well as maximum system pressures with respect to its specifications. The according maximum pressure applies to all modes of operation and cannot be exceeded.

PRESSURE AUTOZERO: You can make an offset correction for the pressure display. Open the purge valve by turning the de-aeration screw, see {4.3} in Fig. 5 on page 7, one turn to the left. Select menu with ➤ and start autozero routine by pressing ▲ or ▼. Offset correction is done automatically.

FLOW: You can enter a correction constant for the selected flow rate in a preset range [1000 - 1400] depending on the type of pump head. This correction changes the number of pump cycles at a selected flow rate. This option allows for balancing different compressibility of solvents.

Only trained users should change that setting!

Select menu with ➤, position cursor on desired numeric field and adjust values with ▲ or ▼.

You have no access to the correction parameter, see Fig. 10. These can be accessed by authorized maintenance personnel only.

PRESSURE LIMITS: Allows to set an individual minimum and maximum system pressure. The pump will switch off automatically if the system pressure exceeds limits. If minimum pressure is set to 0, no check of the minimum pressure is carried out. The absolute maximum pressure is limited by the pump head chosen and by the flow rate selected for this pump head. The detailed pressure limits are given in „Maximum Pressure“ on page 19.

The specified maximum pressure applies to all modes of operation and cannot be exceeded. Units: 0.1 MPa (= bar).

START INPUT: If the START INPUT is set to ENABLED, a program can be started by an external digital signal from the START IN input on the rear side of the pump, see chapter „Connecting the Analog/error/start Socket“ on page 28. The setting STOP Pump will cause the pump to stop pumping, if an external

digital signal is applied on the START IN input. The setting START Pump enables to start pumping by using an external signal.

If the START INPUT is set to START Pump or STOP Pump, no program runs can be started and network control is inactivated.

DATE/TIME: You can use this menu for setting the system date and time of your Preparative Pump 1800.

CONTROL: The instrument can be set to different external control modes. Options are **KNAUER-NET**, where the baud rate is currently limited to 9600, **RS232**, **RS485** and **ANALOG**.

The option KNAUER-NET is used when operating the pump in a network with Knauer instruments and software, see „Using the RS232 Serial Interfaces in the KNAUER-Net“ on page 30. RS232 and RS485 allow for external digital control by a suited terminal program, ANALOG enables external control via the analog inputs, see „Analog control of the pump“ on page 30.

ANALOG OUTPUT: On activation of PRESSURE, the output ANALOG OUT delivers the actually measured pressure value for external monitoring or controlling. Activating SLAVE B, an analog signal is given out via the output that can be used for controlling any second pump, e.g. in a binary high pressure gradient. In this case the control signal (voltage) corresponds to the percentage amount % B , whereas % A will be delivered automatically by the Preparative Pump 1800 itself.

Generally, the voltage range can be chosen between 0-1 V, 0-2 V, 0-5 V and 0-10 V. In any case, the maximum value corresponds to 100 % B (SLAVE B) or to 40 MPa (PRESSURE).

Using the VIEW Menu

The VIEW menu lists existing programs and program links.

An example is given in Fig. 12: Programs 1, 2, 3, 4, 11 and link 21 have been created.


Fig. 12 The VIEW-menu with programs 1-4, 11 and link 21

Using the GLP Menu

The GLP menu reports statistical data for vital parameters of the Preparative Pump 1800. Fig. 13 gives an overview of the GLP menu.


Fig. 13 Structure of the GLP-menu

For an overview of error codes refer to chapter „List of Display Messages“ on page 38.

Using the LINK Menu

The LINK menu reports status information of the currently active program links. It is activated only if a link is running.


Fig. 14 LINK-menu without running link

If a link is active, the menu shows the following structure:


Fig. 15 LINK-menu with running link

In case a link is running, the first line of the menu gives the current status of the links. The field R-xx shows the number of program runs still to be performed. The value of R-xx decreases according to the number of remaining runs.

The second line shows the next step of the link to be performed. Up and down arrows at the right hand side of the display indicate existence and position of additional program steps.

Operating the Preparative Pump 1800

You need to setup your pump prior to operation. Use the information in sections „Menu Structure of the Preparative Pump 1800“ on page 9, „Using the SETUP Menu“ on page 12,

Working Example for setting parameters

For your convenience, this example points out how to configure a binary low pressure gradient using valves A and B.

This example may not comply with your needs and is given for convenience only. You need to select parameters that meet your requirements. Using this example without proper adjustment of parameters may lead to damage of your equipment.

Enter the SETUP menu, see Fig. 10 on page 11.

SOP 3 Selection of correct pump head:

This SOP applies to the Preparative Pump 1800 Firmware Revision 2.xx.

1. Select menu PUMP HEAD by pressing ∇
2. Use \leftarrow or \rightarrow to scroll through the list of pump heads until the installed type of pump appears, e.g. 250 ml, steel. Selection of a pump head automatically adjusts proper absolute maximum system pressure.
3. Press \leftarrow to confirm selection and return to the rhombus.

SOP 4 Selection of operation mode

This SOP applies to the Preparative Pump 1800 Firmware Revision 2.xx.

1. Return to page GRADIENT MODE by pressing \blacktriangle .
2. Press \rightarrow and then \blacktriangle to scroll through the gradient mode options until the wanted appears, e.g. **LPG \rightarrow 220**.
3. Use \rightarrow to position cursor in the second line.

The number of ON's selected refers to all programs in the pumps memory.

4. At any position you can use \blacktriangle or \blacktriangledown to select ON or -- (off).
5. Press the \leftarrow or \rightarrow key as often as you need to position the cursor at the rhombus.
6. Leave this menu page with \leftarrow to return to the main menu.
7. Restart the pump to read out the valve values corresponding to the selected LPG mode .


Fig. 16

Example for a gradient setup with solvents A and B in all programs

If after leaving this item over the rhombus the message **Conflict with existing programs** is displayed, this means that programs with different valve or pump configuration already exist. In this case, either delete these programs according to SOP 16 "Modifying and Deleting Programs" on page 24 or change your valve or pump configuration correspondingly.

Using the Low Pressure Gradient Mode

Two low pressure gradient units (binary A4094 and ternary A4099) are designed for operating with the Preparative Pump 1800. A firmware version of 2.1 or higher is required for the pump.

The valve blocks can be controlled directly in the stand alone mode and with the control software EuroChrom[®], preparative version (release 3.01 and higher).

The ternary LPG valve unit can be used with a flow rate range of 10 up to 220 ml/min. If it is operated with a 500 ml or 1000 ml pump head, it will not operate properly with flow rates above 220 ml/min.

The binary LPG valve unit can be used with a flow rate range of 10 up to 1000 ml/min, mainly in the flow range of 100 up to 750 ml/min.

SOP 5 Installing the ternary LPG valve block

This SOP is valid for the Preparative Pump 1800 Firmware Revision 2.xx and the ternary LPG valve block A4099.

1. Select the **LPG→220** gradient mode in the setup menu and define the channels in use according to SOP 4.

Further installation must be done when the pump is turned off and disconnected from the network!

Operating the ternary valve block with parameter settings of the binary one will destroy the block. Take care, that the pump has the right settings! (Fig. 11 on page 12) If necessary, the pump settings have to be adjusted to the ternary valve block.

2. Remove the eluent inlets from the inlet bushings on the pump head.
3. Install the valve block as illustrated in Fig. 17 on the pump head.
4. In older pump heads leakage's may occur. In this case, insert additional sealings on both bushings.

Make sure that they are inserted above and below the valve blocks mounting base.

5. The eluent suction tubes are connected to the Dynaseal bushings and the three inlets (A, B, and C). If the unit shall be used for a binary gradient, connect a blind fitting to inlet C.
6. Insert the control cable into the LPG valve socket on the back panel of the pump.


Fig. 17 Installing the ternary LPG valve block

SOP 6 Installing the binary LPG valve block

This SOP is valid for the Preparative Pump 1800 Firmware Revision 2.xx and the binary LPG valve block A4094.

1. Select the **LPG→1000** gradient mode in the setup menu and define the channels in use according to SOP 4.

Further installation must be done when the pump is turned off and disconnected from the network!

Operating the binary valve block with parameter settings of the ternary one will lead to malfunction of the valves. Take care, that the pump has the right settings! (Fig. 11 on page 12) If necessary, the pump settings have to be adjusted to the binary valve block.

2. Proceed as described in items 2 to 4 in SOP 5
3. The eluent suction tubes are directly connected to the inlets at the block.
4. Insert the control cable into the LPG valve socket on the back panel of the pump.

Fig. 18 Binary LPG valve block mounted to the Preparative Pump 1800

Using the High Pressure Gradient Mode

If you want to control up to three additional HPLC Pumps K–1800 for a high pressure gradient, you can connect the appropriate number of pumps using the serial RS232 interfaces according to „Using the RS232 Serial Interfaces in the KNAUER-Net“ on page 30.

Activate „GRADIENT MODE“ in the SETUP menu, see „Using the SETUP Menu“ on page 12 and select individual, unique names from HPG A to HPG D for all pumps connected.

The number of pumps in your HPG system can be selected by using the ON option in the second line. Only the pumps selected with ON will be used for calculating the percentage composition in a program.

Every pump automatically takes control over all pumps connected, if execution of a program is started on that pump.

Direct Control of the Pump

If you want to use the Direct Control Mode you have to specify:

- Maximum system pressure
- Flow Rate and
- Solvent Composition if working in LPG or HPG gradient mode

Maximum Pressure

The range for the pressure limits depends on the flow rate selected for the respective pump head. In the lower flow rate range, higher maximum pressures are permitted.

If the flow rate is raised, the maximum pressure limit in the SETUP is adjusted automatically.

Table 2 Maximum pressure limits of different pump heads

Pump Head [ml]	Flow Range [ml/min]	Maximum Pressure [0,1 MPa]
100	0 - 100	400
250	0 - 100	225
	100 - 150	linear reduction 225 → 200
	150 - 250	200
500	0 - 500	100
1000	0 - 350	75
	350 - 600	linear reduction 75 → 50
	600 -1000	50

SOP 7 Setting the Maximum Pressure

In addition to the maximum pressure specified in the SETUP menu, either by the pump head or by the individually set pressure limit, you can specify a different, lower system pressure for the currently active program.

This SOP applies to the Preparative Pump 1800 Firmware Revision 2.xx.

1. Move the cursor to the parameter PRESSURE, see {1.4} in Fig. 2 on page 6.
2. Type in the desired value for the maximum pressure. The selected value is marked with two vertical lines which indicate that this is a programmed, not a measured value.
3. Press any arrow key for confirmation.

SOP 8 Setting the Flow Rate

This SOP applies to the Preparative Pump 1800 Firmware Revision 2.xx. Use this SOP for setting the correct flow rate.

1. Position cursor on field FLOW, see {1.3} in Fig. 2 on page 6.
2. Enter the desired flow rate with the numeric keys.
3. Press any arrow key for confirmation.

SOP 9 Setting the Solvent Composition

This SOP applies to the Preparative Pump 1800 Firmware Revision 2.xx, when working in any gradient mode.

1. Move the cursor to the field % A, see {1.8} in Fig. 2 on page 6.
2. Type in the desired percentage of eluent A. The value in the right hand adjacent field next to your input complements the value to give a total of 100%.
3. Press any arrow key for confirmation.

Events

Events are defined as states or changes of states of the electrical contacts of the Events socket of the instrument. You can either open or close any of the contacts or use a 500 ms pulse. Some times the manual control may be useful, e.g. for test purposes. For connection of cables use sections „Electrical Connections“ on page 28

SOP 10 Setting Events

This SOP applies to the Preparative Pump 1800 Firmware Revision 2.xx.

1. Position the cursor at the event field.
2. Start entry of events by pressing 0, 1, 2 or 3 for the first event. After every entry, the underlined cursor jumps one position to the right. The following table gives the meaning of the possible entries:
 - 0 — Off
 - 1 — On
 - 2 — Pulse
 - 3 — No change with respect to previous status
3. Finish entry of events by pressing any arrow key.


Fig. 19 Relations between event programming and event terminal strip

SOP 11 Start and Stop of Solvent Delivery

This SOP applies to the Preparative Pump 1800 Firmware Revision 2.xx.

1. Apply SOP 8 „Setting the Flow Rate“ on page 19.
2. Use the START/STOP Key, see Pos. {2.1} in Fig. 2 on page 6, to start and stop solvent delivery.

SOP 12 Purging the Preparative Pump 1800

This SOP applies to the Preparative Pump 1800 Firmware Revision 2.xx.

Open the „De-aeration screw for opening purge valve“, Pos. {4.3} prior to starting the purging function to avoid damaging your column system.

1. Purging without opening the purging outlet may cause damage to your column and system.
2. Press purge key {2.2}, see Fig. 2 on page 6, for start and stop of purging. In LPG mode, the pump is purged with the currently active solvent composition.
3. Select a flow rate for purging, see SOP 8 „Setting the Flow Rate“ on page 19.
4. In LPG modes, purging can be performed either in mode MIX, using the currently active solvent composition or separately for each channel A to D. Select the desired mode by positioning the cursor on the appropriate positions MIX, A, B, C or D.

Automatic Piston Backflushing

The Preparative Pump 1800 is equipped with automatic piston backflushing. This allows a longer lifetime of the seals and pistons and removes traces of salt and other decontaminants from the backside of the pistons.

Automatic piston backflushing is effected over a rinsing pump which is installed inside of the Preparative Pump 1800.

SOP 13 Preparing your Preparative Pump 1800 for automatic piston backflushing

This SOP applies to the Preparative Pump 1800 Firmware Version 1.0X or higher. please connect the eluent tubing as follows:

1. Connect the silicone tubing from the accessory kit to position 1 and put the other end into the reservoir bottle (see Fig. 20 below).


Fig. 20 Piston backflushing capillary connections

2. Connect the silicone tubing from the accessory kit to position 4 and put the other end into the waste bottle.

The self-priming piston backflushing pump is equipped with a time control device. Once the Smartline Pump 1000 is powered on, the rinsing pump turns on automatically for 15 seconds and fills the whole tubing and the backside of the pistons with washing solution. The rinsing pump turns on automatically every 2 hours for 15 seconds as long as the Smartline Pump is on. The washing solution flows from the reservoir bottle to the rinsing pump over position 1 and flows out from position 4 through the pump head (positions 2 and 3) and then goes to the waste bottle.

For washing solution we recommend for example a solution of water/methanol (80/20) or pure water.

Please use only pure water, methanol or isopropanol as the washing solution in the rinsing pump. Use of other solvents (for example THF or toluene) can cause damage to the tubing and rinsing pump.

Programming the Preparative Pump 1800

Programs can perform the following actions:

- Determining and setting flow rates
- Determining and setting solvent compositions
- Control of signal events

The Preparative Pump 1800 can store up to twenty programs with a maximum of 100 program lines over all. You can define up to nine links between existing programs.

Features of Programs

Units of the time display are minutes with a decimal representation of seconds, i. e. 0.3 min = 18 s. All programs are saved in the pump's memory and are permanently available. You can enter values with an accuracy of 1/10 min. Time is displayed in 1/100 min steps during program execution.

Programming the solvent composition does not depend on the currently active gradient mode, i. e. HPG or LPG.

New program lines are appended in the last line with the cursor on an asterisk. Program lines can be rearranged manually by changing the time values.

If you generate a new line, all parameters are set to status No Change, displayed as underscore characters, see Fig. 22 on page 25, part D. The No Change-feature allows changing single parameters without reentering the values for others.

When parameters have different values in different time lines, their value will be calculated and changed linearly in the time span in between. This is important e.g. if programming linear gradients.

SOP 14 Creating a Program

Use this SOP for entering program information in your Preparative Pump 1800. Use „Features of Programs“ above.

1. Select the operation screen (main menu).
2. Position cursor on field {1.2} EDIT PROG
3. Enter number of desired program.
4. Press any arrow key to enter edit mode for new program.

The cursor is positioned at Time 000.0 (fixed start time). Trying to change the start time will produce the message Start Time Is Fixed in the display.

5. Position cursor on desired field.
6. Enter desired values in the desired fields according to SOP 8, „Setting the Flow Rate“, SOP 9 „Setting the Solvent Composition“ and SOP 10 „Setting Events“ on page 20.
7. Position cursor on time field and press u for creation of an additional time step.

SOP 15 Execution of Programs

This SOP applies to the Preparative Pump 1800 Firmware Revision 2.xx.

1. Load the program you want to execute by entering the program number in the LOAD PROG-field. Confirm with any arrow key.
2. Press START key to start the pump. The display changes to **ON** and the pump begins continuous operation with the parameters specified in the first program line.
3. Position cursor on field **ON**
4. Start the program by pressing the numeric key 1. The display changes to **RUN**.
5. During execution of programs you can press key 0 to freeze the current status of the pump. The display changes to **HOLD** and the field TIME is blinking.
6. If you press 1 you can activate program execution again or you can press **STOP** to abort execution.

During program execution the current solvent composition can be observed in fields % A and % B. At the end of the program the pump returns continuous pumping with the parameters of the first line, ON is displayed.

SOP 16 Modifying and Deleting Programs

This SOP applies to the Preparative Pump 1800 Firmware Revision 2.xx. Use this SOP for modifying and deleting programs.

1. Enter the number of the program to be modified in field EDIT PROG.
2. Move the cursor to the value you want to modify.
3. Enter new values. Old ones will be overwritten. Wrong entries can be corrected the same way.
4. Confirm your entries by pressing any arrow key.

Delete program lines as follows:

5. Enter 0 (zero) in the TIME FIELD and press an arrow key. The message **Delete this line? Confirm by cursor** appears. Pressing any arrow key deletes the program line. Pressing any numeric key or none at all for several seconds aborts the process of deletion.

Deletion of programs or program links:

6. Select field EDIT PROG.
7. Enter 0 in the field EDIT PROG
8. Confirm deletion according to step 5 of this SOP.

SOP 17 Using a Wake-Up Program

This SOP applies to the Preparative Pump 1800 Firmware Revision 2.xx.

The Preparative Pump 1800 can execute a program at a preset time.

Make sure that date and time of the Preparative Pump 1800 are set up properly according to section „Using the SETUP Menu“ on page 13.

A Wake-Up Program is created as follows:

1. Enter the reserved program number 30 in the field EDIT PROG.
2. Enter the desired program for Wake-Up in field P: and the desired time and date in fields at: and on:.
3. Load reserved program 30.
4. If necessary, switch pump to OFF with the START/STOP key.

While waiting for the programmed time, the display is blinking. The pump will start program execution with the selected program at the preset time and switch into RUN mode.

5. Execution of Wake-Up Program can be aborted by entering the number of another program.


Fig. 21 Display of Wake-Up program

Example for creating a program

Generation of program no. 7 with the following functions: Constant flow rate of 100.0 mL/min. Linear gradient changing composition from 90 % A and 10 % B at start time to 10 % A and 90 % B at time 20 minutes.

This program is based upon the configuration in section „Operating the Preparative Pump 1800“ on page 15. This setup must be performed before program definition. Remember, underscored characters indicate that the corresponding parameter is continuously calculated.

1. Select the operation screen
2. Move cursor to the field Edit Prog, press 7 to generate program no. 7.
3. Press ∇ to enter the Edit mode for the new program. The cursor is positioned at Time 000.0.
4. Press \blacktriangle to set the flow rate: Press 1, 0, 0 and then ∇ to enter the value 100.0 mL/min. The display should look like part A see Fig. 22, cursor blinking in the TIME field.


Fig. 22 Display of single steps, when creating a program

5. Press \blacktriangleright to enter field % A , then 9 and 0 and \blacktriangleleft to confirm the input. The amount for % B is automatically calculated and indicated as 10 %. The display should now look like part B in the figure.
6. Press now ∇ again to create the next program line by defining a new time step. The display now should look like part C in the figure, the cursor blinking on the first asterisk.
7. Now press 2 and 0 and then \blacktriangleright to define time step 20.0 min as well as to enter field % A.
8. Press 1 and 0 for 10 % A and then \blacktriangleleft to confirm the input. The amount for % B is automatically calculated and indicated as 90 %. Part D in the figure shows the current situation.

After loading, the program can now be executed.

Using Program Links

Program links are combinations of existing programs, which can be created and edited like any programs.

Links use reserved program numbers **21 ... 29**. These numbers are automatically converted like **21 = L1, 22 = L2** etc. The Preparative Pump 1800 can store up to 9 different links.

Links use line numbers from **No01 to No99**, each line referencing a linked program. The parameter **R** demands the number of runs to be performed with that program. A wait status can be programmed. The maximum number of programs referenced in all links is limited to 100.

Any program can be referenced multiple times in a link. A program link can not reference other links.

SOP 18 Creating Program Links

This SOP applies to the Preparative Pump 1800 Firmware Revision 2.xx. Use this SOP for creating program links in your Preparative Pump 1800.

You can create links between existing programs as follows:

1. Enter desired number for link in field EDIT PROG and press **➤** for confirmation.
2. Press any numeric key to create a line number. Lines are numbered automatically and can be interchanged by renumbering.
3. Use **➤** to move the cursor to field **P** and choose the desired program number.
4. Press **➤** again to select field **R**.
5. Enter number of runs to be performed for the selected program.
6. Select field **Wait** and enter values 1 or 0.

If Wait is set to 1, the pump waits for an external start signal or an user pressing numeric key 1 before this line is executed. On waiting the pump switches to wait-status.

If Wait is set to 0, the program lines of the link are processed continuously.

Example for a program link

Part A in Fig. 23 shows the first step, after entry of link number, cursor blinking at **No****. Part B in the figure shows the programming of the first line.


Fig. 23 Display while editing a program link

During execution of a link the LINK menu is available in addition to the common display. In the LINK menu the momentary status of the running link is shown. Move cursor to the LOAD PROG-field and press < to enter this menu page. For more information refer to Fig. 10 on page 17.

Press > to exit link menu.

SOP 19 Executing Program Links

This SOP applies to the Preparative Pump 1800 Firmware Revision 2.xx.

1. Enter the desired link number into the field LOAD PROG and confirm by pressing any arrow key.
2. Press START to start the pump. The display switches to ON.
3. Position cursor on field ON.
4. Start the link by pressing numeric key 1. The display switches to RUN.

After execution of a link the Preparative Pump 1800 stops solvent delivery, the display switches to OFF.

Connecting the Preparative Pump 1800 to other Instruments

Hydraulic Connections

Make sure that all hydraulic connections are suited for the system pressure and flow rate of your HPLC System.

100, 250 and 500 ml Head

Solvent tubings are connected to the pump head according to Fig. 24. Make sure that the flat side of the ferrule is directed to the pump head and that the fitting screw is fastened hand tight. According to the local conditions the tubings can be connected either to the inlets valve housing (A) or to the eluent inlet (B).


Fig. 24 100, 250 and 500 ml pump head: hydraulic connections

1000 ml Head

Push the tubes directly over the connectors, and secure them by means of appropriate clamps.

Fig. 25 1000 ml pump head: hydraulic connections

The pump should be rinsed with a methanol-water mixture (1 : 9 v/v) before first use. Make sure that all connections are suitable for the system pressure and flow rate of your liquid transportation system.

Electrical Connections

Using the EVENT and ANALOG/ERROR/START Sockets

Electrical signal lines from other instruments can be connected to

- the Analog/error/start or
- Events socket or
- RS232 serial interface

on the rear side of the Preparative Pump 1800, see positions {5.2},{5.4}, and {5.5} in Fig. 6 on page 8.

Events can be controlled in the main menu or in programs. You can use SOP 10 „Setting Events“ on page 20.

Do not touch the electrical contacts of the socket lines. Electrostatic discharges can lead to damage of the electronics of the Preparative Pump 1800.

Configuring Plug Strips

For external control of the Preparative Pump 1800 you have to configure a special plug strip according to the local conditions. You can use the plug strips with 4 or 10 connectors, enclosed in the accessories, for easy installation of electrical connections.

SOP 20 Mounting Plug Strips

Use this SOP for connecting wires to the plug strips.


Fig. 26 Mounting Plug Strips

Connecting the EVENTS socket

You can control the states of individual contacts of the EVENTS socket using SOP 10 „Setting Events“ on page 20 an.

Specifications of the Event signals

TTL (Transistor-Transistor-Logic) connections: TTL compatible outputs (max. consumption: 40 mA, $U_{CEmax} = 24\text{ V}$) referred to GROUND.

Table 3 Technical specifications of the Event signals

State	Control outputs	OC (Open Collector)-connections	Relay connections
0	Low (< 500 mV)	Inactive	Contact open
1	High (> 2 V)	Active	Contact closed
2	Pulse, > 2 V for 500 ms	Pulse, active for 500 ms	Pulse: Contacts 500 ms closed

Connecting the Analog/error/start Socket

The analog/error/start terminal strip serves to receive or send start or error signals from/to other instruments. Additionally, a further voltage signal is available for monitoring the system pressure or for controlling a slave pump, see item ANALOG OUTPUT in „Using the SETUP Menu“ on page 12. The ANALOG IN input allows external analog control, see „Analog control of the pump“ on page 30.


Fig. 27 The REMOTE terminal strip

Table 4 Technical specifications of analog/error/start signal connections and picture of connector

Signal	Description
Start connections	
START IN	Is activated by a 0 Volt or low signal or a short circuit $I_{IN} = 10\text{ mA}$
START OUT	OC (open collector) output is active for 500 ms
Error- connections	
ERROR IN	When receiving a 0 V (error) signal, the message Error signal was detected is displayed and the pump is stopped; $I_{IN} = 10\text{ mA}$
ERROR OUT	OC output — active until the error condition is removed
Analog- connections	
ANALOG OUT	Analog voltage signal — Gives either measured system pressure (1 V per 40 MPa; Offset for 0 MPa is some mV) or a control voltage for slave pump B. Adjustable voltage range between 1, 2, 5, and 10V maximum output.
ANALOG IN	Analog input signal controls flow rate; default 1 V per 10 ml/min for 100 ml head, 1V per 25 ml/min for 250ml head, 1V per 50 ml/min for 500ml head, and 1 V per 100ml/min for 1000 ml head. Maximum input voltage is 10 V.

Analog control of the pump

When the option ANALOG is activated, see item CONTROL in „Using the SETUP Menu“ on page 12, the Preparative Pump 1800 can be controlled by an external voltage signal.

After selecting ANALOG, an offset correction for the input voltage can be set. Move cursor to Off's and press ∇ or \blacktriangle for automatic offset correction of the input voltage.

In the second line of the Analog Control window the current flow rate is indicated. The flow rate is calculated according to the voltage applied at the input. The ratio between input voltage and flow rate can be varied as follows:

Apply a voltage of 10V maximum to the connectors **Analog Ground** and **Analog In**. The corresponding flow rate will be displayed. Move the cursor to that field. Enter the flow rate desired for the given voltage. Obey the limitations given by the capacity of the pump head and the maximum voltage of 10 V. Confirm the input pressing the ∇ or \blacktriangle button. Entering 0 (zero) reactivates the default settings. The allowed maximum flow rates for 1 V are given in the following table:

Pump Head	Flow _{Max} for 1V Input Voltage
100ml	15,0ml
250ml	37,5ml
500ml	75,0ml
1000ml	150,0ml

In case of these settings the maximum input voltage is for all pump heads 6.66 V. Similarly the maximum flow rates will be reduced, if smaller flow rates are defined for 1 V, e.g. setting the flow rate for 1 V on 1 ml for a 100 ml pump head results in a maximum flow rate of 10ml/min at 10 V.

Using the RS232 Serial Interfaces in the KNAUER-Net

This chapter will explain the necessary connections for operating the Preparative Pump 1800 controlled by EuroChrom[®] or ChromGate[®]. For detailed information refer to the user manual of the software.

The two RS232 serial interfaces on the rear side of the device, see {5.4} in Fig. 6 on page 8, enable digital data transfer between the Preparative Pump 1800 and other instruments like other pumps, detectors or a PC which runs a software package like ChromGate[®] or EuroChrom[®], Preparative Version. All the devices are communicating with each other and form an integrated network. In the Setup Menu control by KNAUER NET has to be chosen.

All devices are connected in a token ring in the following way: The first RS232 interface (the left one) of the first device is connected to the second RS232 interface (the right one) of the next device and so on. The ring is closed by connecting the RS232 interface of the last device to the free interface of the first one, see Fig. 28.

Fig. 28 Connection scheme of three Preparative Pumps 1800

For connections with a PC you need a Y-shaped cable which is part of the EuroChrom 2000 Software/Hardware package. You may include also other network integratable KNAUER instruments, like the fast scanning Spectro Photometer K-2600 or the KNAUER interface box.


Fig. 29 Connection scheme of a Preparative Pump 1800 with PC and other instruments

Easy Maintenance

Changing and Cleaning the Pump Head

SOP 21 Removing the Pump Head

This SOP refers to the Preparative Pump 1800.

For exchange of the pump head, or for disassembling it in order to clean valves or replace piston seals, piston rods etc., the pump head can easily be removed.

1. Purge the pump head with a suitable cleaning agent and then with distilled water. You can use SOP 12 „Purging the Preparative Pump 1800“ on page 21.

If organic solvents remain in the pump head, danger of skin irritation may exist.

2. Remove the two solvent tubings from the „Eluent inlets“, Pos. {3.2} see Fig. 5 on page 7.
3. Loosen the two „Eluent outlets to transducer“, Pos. {3.7} and the two „Inlets to transducer“, Pos. {4.1}.
4. Remove the two „Connection capillaries“, Pos. {3.8}, beginning with the upper „n“-shaped one.
5. Loosen the „Inlets valve housing“, Pos. {3.5} as well as the fittings at „Outlets valve housing“, Pos. {3.6}.
6. Using a hexagonal spanner no. 4 (4 mm), loosen just two diagonally opposed „Head set screws“, Pos. {3.7}. Remove the screws. (see Fig. 30 on page 33).
7. Carefully loosen the two remaining screws, alternating from one to the other, approx. half a turn. This prevents the pump head from tilting and becoming damaged. Once the spring tension has been reduced, hold the pump head firmly in one hand while removing the screws completely with the other hand.
8. Carefully remove the pump head.

SOP 22 Removing and Checking Piston Rods

This SOP refers to the following pump heads:

- 100 ml stainless steel, KNAUER order number A 4029(-1)
- 250 ml stainless steel, KNAUER order number A 4021(-1)
- 500 ml stainless steel, KNAUER order number A 4038(-1)
- 1000 ml stainless steel, KNAUER order number A 4022(-1)

1. Remove Pump head as described in SOP „Removing the Pump Head“.

If you only wish to check the piston rods, you do not need to disassemble the pump head any further.

2. The „Piston rods“, Pos. {3.17} in Fig. 30 on page 33 may be removed using pliers. Take the tip of the piston using the pliers, and pull it out carefully in a straight line.

If the rods are broken, you must check the entire pump head for damage.

SOP 23 Disassembling the Pump Head I

This SOP refers to the pump heads 100 ml, 250 ml, 500 ml, and 1000 ml in compact design.

All positions of components refer to Fig. 30 „Explosion diagram of pump heads, compact design“ on page 33.

1. Apply SOP 22 "Removing and Checking Piston Rods" on page 32.
2. Loosen the two „Retaining plate screws“, Pos. {3.19} half a turn, alternating from one to the other to avoid damaging the „Retaining plate“.
3. Because the two screws are very tight, it may be helpful to either clamp the pump head or to press one of its side surfaces against a table with one hand while loosening the screws.
4. Unscrew the two screws strictly alternating due to strong force of the „Springs“, Pos. {3.15} behind the plate, and remove them.
5. Remove the „Retaining plate“, Pos. {3.18}.
6. You can now remove the „Guides for springs“, Pos. {3.16}, the „Springs“, Pos. {3.15} as well as the two „Distance plates“, Pos. {3.13}.
7. Using the SW 10 spanner, loosen and remove the „Spacing bolts“, Pos. {3.14}. These bolts are seated very tightly. Follow the tip given in step 3.
8. Lift and remove the „Pump head housing, basic plate“, Pos. {3.B}.
9. To remove the piston seals, Pos. {3.11} in the basic plate {3.B} and in the front block {3.A}, use an appropriate tool, e.g. a screw driver and carefully lever the four piston seals.


Fig. 30 Explosion diagram of pump heads, compact design


Fig. 31 Explosion diagram of pump heads, inlay design

3.A Pump head housing, front block	3.10 Check valves, inlet and outlet
3.B Pump head housing, basic plate	3.11 Piston seals
3.1 Piston backflushing capillaries 1/16"	3.12 O-shaped sealings
3.2 Eluent inlets	3.13 Distance plates
3.3 Seal holder	3.14 Spacing bolts
3.4 Pump head inlays	3.15 Springs
3.5 Inlet valve housing	3.16 Guides for springs
3.6 Outlet valve housing	3.17 Piston rods
3.7 Pump head set screws	3.18 Retaining plate
	3.19 Retaining plate screws

SOP 24 Disassembling the Pump Head II

This SOP refers to the pump heads 100 ml, 250 ml, 500 ml, and 1000 ml in inlay design.

All positions of components refer to Fig. 31 „Explosion diagram of pump heads, inlay design“ on page 34.

1. Perform steps 1 through 5 of SOP 23 „Disassembling the Pump Head I“.
2. Remove the outer piston seals {3.11} through the basic plate {3.B} using the piston from the inlay {3.4}.
3. Loosen the spacing bolts {3.14} with a spanner SW 10. Because the two screws are very tight, it may be helpful to either clamp the pump head or to press one of its side surfaces against a table with one hand while loosening the screws.
4. Remove the pump head basing plate {3.B}.
5. Now the inlays {3.4} can be easily removed from basic plate.
6. To remove the inner piston seals, Pos. {3.11} in the seal holder {3.3}, use an appropriate tool, e.g. a screw driver and carefully lever the piston seals.

7. To remove the seal holders {3.3} from the pump head front housing {3.A} first remove the inlet and outlet valve housings {3.5 and 3.6}.

SOP 25 Assembling the Pump Head and Exchanging Piston Seals I

This SOP refers to the pump heads 100 ml, 250 ml, 500 ml, and 1000 ml in compact design.

All positions of components refer to Fig. 30 „Explosion diagram of pump heads, compact design“ on page 33.

1. It is recommended to exchange all four piston seals located in the pump heads after disassembling the pump head.
2. With the open side facing downwards, carefully press the „Piston seals“, Pos. {3.11} into the „Pump head housing, front block“, Pos. {3.A} and „Pump head housing, basic plate“, Pos. {3.B}, making sure to keep them straight.

Fig. 32 Open side of piston seal

3. Insert the O-ring sealings part {3.12}.
4. Put the two housing parts together, positioning the O-ring side of part {3.B} against part {3.A}. Be sure that the backflushing capillaries are positioned on the top side of the head - see label on front of part {3.A} for orientation.
5. Insert the two „Spacing bolts“, Pos. {3.14}, screw them in with the hand and tighten them as firmly as before using the SW 10 spanner.
6. Put in the two „Distance plates“, Pos. {3.13}, position the „Springs“, Pos. {3.15} and the „Guides for springs“, Pos. {3.16}.
7. Put the „Retaining plate“, Pos. {3.18} over the „Guides for springs“, Pos. {3.16}, with the conical bore of the plate directed outwards.
8. Insert the two „Retaining plate screws“, Pos. {3.19}, press the plate down with one hand and tighten the screws with a screw driver, strictly alternating due to strong force of the „Springs“, Pos. {3.15}.
9. Carefully insert the piston rods, see „Piston rods“, Pos. {3.17} without bending or quenching the rods.
10. The „Retaining plate screws“, Pos. {3.19} must be tightened that they are seated as securely as before.

SOP 26 Assembling the Pump Head and Exchanging Piston Seals II

This SOP refers to the pump heads 100 ml, 250 ml, 500 ml, and 1000 ml in inlay design.

All positions of components refer to Fig. 31 „Explosion diagram of pump heads, inlay design“ on page 34.

1. It is recommended to exchange all four piston seals located in the pump heads after disassembling the pump head.
2. With the open side facing downwards, carefully press the „Piston seals“, Pos. {3.11} into the seal holders {3.3} and pump head inlays {3.4}, making sure to keep them straight.

3. Insert the seal holders {3.3} into the pump head front housing {3.A} and tighten them by mounting the the inlet and outlet valve housings {3.5 and 3.6}. The inlets have to be placed below the inlay label and the outlets above of them.
4. Insert the pump head inlays {3.4} into the pump head basic plate. Occasionally exchange the O-rings {3.12}.
5. Put he second pair of O-rings into the groove of inlay and basic plate.
6. Assemble the housing parts with outlets and back flushing capillaries upwards.
7. Insert the two spacing bolts {3.14}, screw them in with the hand and tighten them as firmly as before using the SW 10 spanner.
8. Put in the two „Distance plates“, Pos. {3.13}, position the „Springs“, Pos. {3.15} and the „Guides for springs“, Pos. {3.16}.
9. Put the „Retaining plate“, Pos. {3.18} over the „Guides for springs“, Pos. {3.16}, with the conical bore of the plate directed outwards.
10. Insert the two „Retaining plate screws“, Pos. {3.19}, press the plate down with one hand and tighten the screws with a screw driver, strictly alternating due to strong force of the „Springs“, Pos. {3.15}.
11. Carefully insert the piston rods, see „Piston rods“, Pos. {3.17} without bending or quenching the rods.
12. The „Retaining plate screws“, Pos. {3.19} must be tightened that they are seated as securely as before.

SOP 27 Installing the Pump Head

This SOP refers to the following pump heads:

- 100 ml stainless steel, KNAUER order number A 4029(-1)
- 250 ml stainless steel, KNAUER order number A 4021(-1)
- 500 ml stainless steel, KNAUER order number A 4038(-1)
- 1000 ml stainless steel, KNAUER order number A 4022(-1)

1. Make sure that the pump head is properly assembled.
2. Position the head in a straight line onto the pump housing.
3. Tighten all four „Pump head set screws“, Pos. {3.7} a few turns by hand.
4. Alternating from one to the next, tighten two diagonally opposed screws half a turn at a time, until the pump head is correctly seated.
5. Tighten the two other screws.
6. Make sure that all four „Pump head set screws“, Pos. {3.7} are securely tightened.
7. Install the two „Connection capillaries“, Pos. {3.8}, beginning with the lower „s“-shaped one, and firmly tighten all fittings.

SOP 28 Cleaning and Replacing Check Valves

This SOP refers to the following pump heads:

- 100 ml stainless steel, KNAUER order number A 4029(-1)
- 250 ml stainless steel, KNAUER order number A 4021(-1)
- 500 ml stainless steel, KNAUER order number A 4038(-1)
- 1000 ml stainless steel, KNAUER order number A 4022(-1)

If the check valves become dirty they will no longer open and close correctly. You can remove the check valves for cleaning from the pump head after loosening and removing the „Inlet valve housing“, Pos. {3.5} and „Outlet valve

housing“, Pos. {3.6}, and then disassemble them according to the following instructions.

1. Place the valves in a suitable cleaning solution. Use an ultrasonic bath to clean the valve. If malfunction persists, use steps 2...5 of this SOP.
2. Using a knife or similar, carefully remove the valve seal from the housing; see Fig. 23 on page 41.
3. Remove the valve ball out of the housing.
4. Clean the individual parts. We recommend an ultrasonic bath.
5. Assemble the check valves in reverse order, use Figure Fig. 33.


Fig. 33 Components of the check valves of all pump heads available

List of Display Messages

The following list shows messages that may appear in the main screen of the Preparative Pump 1800 and the pages on which they are described.

SYSTEMTEST OK	9
Conflict with existing programs	16
Start Time Is Fixed	23
Delete this line? Confirm by cursor	24
Error signal was detected	29

List of Error Codes

Table 5 Error codes in GLP Menu

Error Code	Meaning
5	Motor failure
9	ERROR IN signal detected
10	Switch off due to exceeding maximum current
12	Switch off due to exceeding maximum pressure
13	Switch off due to falling below minimum pressure
23	Motor blocked
33	Calibrating values re-initialized
34	Calibrating curve re-initialized
35	Combination from 33 + 34

Packing List

The Preparative Pump 1800 will be delivered together with one of the following pump heads:

100ml pump head, order number A50501,
250ml pump head, order number A50511,
500ml pump head, order number A50521, or a
1000ml pump head, order number A50531.

The standard version will be delivered with stainless steel inlays. For biocompatible working also pumpheads with titanium inlays are available. The order numbers of pumps equipped with such pump heads bear the extension V1 (with KEL-F seals) or V2 (with PEEK seals).

The delivery consists of::

- Preparative Pump 1800 basic unit with pump head installed
- User Manual
- Power supply cable
- Power cable for connecting further instruments
- RS232 cable
- Set of plug strips and signal cables including
 - 3 plug strips with four connections
 - 1 plug strips with ten connections
 - 3 m signal cable with two poles
 - 1 m signal cable with one pole
- Set of tools including
 - 1 wrench 12 x 17 mm
 - 2 wrench 8 10 mm
 - 2 wrench 1/4'' x 5/16''
 - 1 hexagonal spanners 4 mm
 - 1 cross screw driver no. 2
- Set of fuses 5A for 250 V
- Syringe set including
 - 1 syringe
 - 1 Luer-Lock needle
 - 0,1 m PTFE tube (ID 1.5 mm; OD 2.1 mm)
- Eluent outlet set
 - 2 fitting screws
 - 2 ferrules
- Accessory kit
 - for 100 ml, 250 ml or 500 ml pump head including:
 - 4 tube fitting screws
 - 6 tube ferrules
 - 2 x 1 m PTFE tubing (ID 3 mm; OD 4 mm)
 - 0,5 m PTFE tubing (ID 1.5 mm; OD 2.1 mm)
 - for 1000 ml pump head including:
 - 2 x 1 m PTFE tubing (ID 7 mm; OD 9 mm)
 - 0,5 m PTFE tubing (ID 1.5 mm; OD 2.1 mm)

Spare Parts and Accessories

	KNAUER Order Number
User Manual	V 7116
Set of plug strips and signal cables	A 14022
Optional Pump Heads	
Pump Head 100 ml incl. accessory kit	A 4029-1
Pump Head 250 ml incl. accessory kit	A 4021-1
Pump Head 500 ml incl. accessory kit	A 4038-1
Pump Head 1000 ml incl. accessory kit	A 4022-1
Spare Parts for Pump Head 100 ml see „Explosion Diagram of 100 ml, 250 ml, 500 ml and 1000 ml pump heads“ on page 34	
Check valve {3.10}	A 1122
Piston rod 1/4" {3.17}	A 0747
Set of Piston seals including	
2 piston seals for 1/4" pistons{3.11}	
1 O-shaped sealing {3.12}	A 0746
Set of Hydraulic connections including	
10 tube fitting screws and 10 tube ferrules	A 2003
Inlet junction including	
2 Inlet valve housings	
1 T piece, fitting screw with ferrule	
1 left-hand and 1 right-hand connection tube	A 1121
Spare Parts for Pump Head 250 ml	
Check valve {3.10}	A 1122
Piston rod 3/8" {3.17}	A 1017
Set of Piston seals including	
2 piston seals for 3/8" pistons{3.11}	
1 O-shaped sealing {3.12}	A 1046
Set of Hydraulic connections including	
10 tube fitting screws and 10 tube ferrules	A 2003
Inlet junction including	
2 Inlet valve housings	
1 T piece, fitting screw with ferrule	
1 left-hand and 1 right-hand connection tube	A 1121
Spare Parts for Pump Head 500 ml	
Check valve {3.10}	A 1080
Piston rod 9/16" {3.17}	A 1478
Set of Piston seals including	
2 piston seals for 9/16" pistons{3.11}	
1 O-shaped sealing {3.12}	A 1479
Spare Parts for Pump Head 1000 ml	
Check valve {3.10}	A 1080
Piston rod 3/4" {3.17}	A 1015
Set of Piston seals including	
2 piston seals for 3/4" pistons{3.11}	
1 O-shaped sealing {3.12}	A 1077

Technical Data

Pump Head independent	
Delivery System	Double Piston Pump with Al ₂ O ₃ Pistons
Flow Rate Precision	± 2 % (within 2 - 50 % of flow range; Methanol: Water 80 : 20)
Reproducibility	0.1%
Pulsation	Electronic reduction of residual pulsation
System Protection	Soft Start, P min and P max adjustable
Control	2 RS 232 interfaces, 8 event outputs (TTL, OC, relays), remote connections (start, error), analog pressure monitor, analog input (flow rate)
Storage	Up to 20 programs and 9 program links
Gradient	Expandable to a High (HPG) or a Low Pressure Gradient (LPG) system with up to 4 eluents HPG: Configurable master pump system with up to 3 slave pumps LPG: Control of up to 4 magnetic valves 24 V, 5 W
Power Supply	115/230 V, 47 - 63Hz, 450 W
Dimensions	220 x 185 x 340 mm (W x H x D)
Weight	10.1 kg

Pump Head depending				
Pump Head	100 ml	250 ml	500ml	1000 ml
Flow Rate Range (ml/min)	0.1 - 99.9	0.1 - 249.9	0.1 - 499.9	0.1 - 999.9
Maximum Pressure (bar)	400	200	100	50
Piston Stroke (mm)	8.3	8.3	8.3	8.3
Stroke Volume per Cycle (ml)	0.53	1.18	2.66	4.73
Diameter (mm)	6.35	9.525	14.29	19.05

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SYSTEMTEST OK	48
Conflict with existing programs	55
Start Time Is Fixed	63
Delete this line? Confirm by cursor	64
Error signal was detected	69

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Warranty statement

The warranty period of the Preparative Pump 1800 is 12 months beginning from the date of dispatch from Berlin. Operation inconsistent with manufacturer's instructions or damage caused by unauthorized service personnel are excluded from guarantee. Damage caused by blockages and wear and tear parts such as fuses and seals are not covered by the guarantee. Defective pumps should be sent to the manufacturer for repair.

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If we find a defect covered by the guarantee, repair or replacement, at our discretion, will be carried out free of charge. Packing and transport costs are borne by the purchaser.

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Declaration of conformity

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**Manufacturer's name and address:
Herstellername und -adresse**

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**Preparative Pump 1800
with different pump heads**

Order Numbers, Bestellnummern: **A 505xx**

complies with the following requirements and product specifications:

- Low Voltage Ordinance (73/23/EWG)
EN 61010-1 (08/2002)
- Engineering Guidelines (89/392/EWG)
- EMC Ordinance (89/336/EWG)
EN 6100-3-2 (10/2006)
EN 61326-1 (10/2006)

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INDEX (english)

A			
accessories	40		
analog control	30		
arrow keys	10		
C			
check valve	36		
connection			
additional pumps	31		
electrical	28		
hydraulic	27		
other instruments	31		
control	9		
cursor	9		
D			
data entering	10		
date	13		
direct control	19		
display messages	38		
E			
Error code	38		
event signal	29		
Events	20		
external control	13		
F			
flow range	3		
flow rate			
correction constant	12		
range	41		
setting	19		
G			
gradient mode	12, 16		
H			
HPG	12		
K			
KNAUER-Net	30		
L			
link	15		
LPG	12, 16, 51		
M			
menu			
GLP	14		
link	15		
selection	10		
setup	11		
		structure	9
		view	13
O			
		operation mode	16
P			
		parameter setting	15
		piston backflushing	21
		piston rod	32
		piston seal	35
		plug strip	28
		precision	41
		pressure	
		autozero	12
		limit	12, 19
		program	22
		deleting	24
		editing	23, 25
		execution	23
		features	22
		link	26
		modifying	24
		wake up	24
		pulsation	41
		pump head	
		assembling	35
		connection	7, 27
		disassembling	33, 34
		explosion diagram	33, 34, 35
		identifying	6
		installation	36
		purging	21
		removing	32
		selection	12, 16
R			
		RS232 serial interface	30
S			
		solvent composition	20
		spare part	40
		start and stop	21
		start input	12
		start time	23
T			
		technical data	41
		time	13
V			
		valve	12
		valve block installation	17, 18

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