



**Note:** Read the corresponding technical documentation to use the correct filter.

## Changing filter cartridges in biocompatible pumps with an installed PEEK filter (M2499E)

**Note:** This type of filter releases particles that might lead to fluidic blockages. Therefore, this filter is no longer available.



Replace the PEEK filter with the titanium filter (G1452A) or with the empty cartridge (A9652) together with the black bushing (P3830H) provided with the kit.



**Titanium filter**  
(G1452A + P3830H)



**Empty cartridge**  
(A9652 + P3830H)

Titanium filters (G1452A) are available in a set of three (order number A9661).

