

Overview power plugs and interfaces

1. Allocation of power plug types to devices

Every device is supplied with a power plug of the AZURA® series (cold-device plug) in the suitable country-specific version.

1.1 Overview of country-specific power plugs, routers and switches

If no suitable adapter is available for a specific country, contact the responsible distributor: www.knauer.net/en/Support/Distributors-worldwide

Power plugs/routers/switches	Article no. USA	Article no. UK	Article no. Switzerland	Article no. Europe	Article no. Argentina	Article no. China	Article no. Australia
Power plug AZURA® series (cold-device plug)	M1651	M1278	M1597	M1642	M3233	M3381	M3439
Power plug Smartline series	M1279	M1277	M1479-1	M1479	-	-	-
Router (power plug incl.):	A64809INT	A64809INT	-	A64809	-	-	-
▪ MikroTik	M1651	M1278	M1597	M1642	M3233	M3381	M3439
▪ NetGear, D-Link	Adapter: M0447V2	Adapter: M0447V1					
Switch (power plug incl.):							
▪ Switch 8x LAN	A3119INT Adapter: M0447V2	A3119INT Adapter: M0447V1	-	A3119	-	-	-
▪ Switch 16x LAN	A3129INT Adapter: M0447V2	A3129INT Adapter: M0447V1	-	A3129	-	-	-

 **Note:** To connect multiple devices, we provide a special power plug for up to four AZURA® devices (Europe), article no. A12345.

1.2 Exceptions of allocation

Device	Power plug type
<ul style="list-style-type: none"> BlueShadow Pump 40 P BlueShadow Detector 40D/50D Smartline Degasser (article no. A5328) Osmometer 	Smartline series (see table in chapter 1.1)
<ul style="list-style-type: none"> Router Switch 	Power plug is supplied. For outside Europe, a suitable adapter or cable is supplied (see table in chapter 1.1).
<ul style="list-style-type: none"> Degasser (article no. AZE03, AZE03-1, AZE02-1) 	Power plug is supplied for US, UK, Europe, Australia.
<ul style="list-style-type: none"> Pressure control (article no. AZG10) Pressure sensor (article no. AZG10-1) Airsensor (article no. A70092, A70093, A70082) Interface Box (article no. AZB00XA) 	Power distributor (article no. AZS80SA) and accessories kit with 1x power plug (article no. F1518) is needed. The distributor can provide power for up to 6 devices. Only one power distributor per system is required. Power plug for China or Australia: <ul style="list-style-type: none"> Article no. M3027D (China) Article no. M3027C (Australia)
<ul style="list-style-type: none"> Tablet for Mobile Control 	The tablet includes an european power plug.

2. Allocation interfaces to devices

Currently, desktop PCs from KNAUER have one serial interface (RS-232, DE9). If more than one device with a serial interface has to be connected, install an additional serial interface for the computer. Nearly all laptops have no serial interface. Therefore, in general a serial interface has to be installed.

To operate the following devices, use the USB to serial adapter (article no. A3108) for a single device (also works with Shimadzu RF-20A/Axs), or the USB 4COM (article no. A3114) for more than one device:

- Sedex 85 LT
- Shimadzu RF-20A/Axs
- Bronkhorst Flowmeter
- GJC Flowmeter
- Chiralyser-MP
- GABI* Gamma Spectrometer
- Osmometer (only with software)



Note: If the tablet for Mobile Control should be connected via LAN and not via WLAN, the USB-to-LAN Adapter (article no. A96181) is required. You find the driver on the KNAUER website: www.knauer.net/en/usb-lan-adapter.

Appendix: Detailed overview of devices by power plug type

Power plug AZURA® series (cold-device plug)

All devices of AZURA® series

PCs and monitors

Preparative pumps BlueShadow 80P (article no. APD20xx)

Micro devices

- BlueShadow Pump 10P/20P
- BlueShadow Detector 10D
- Degasser 20DG (article no. AZE02)

Detectors

- RF20A (article no. A59200)
- RF20AXS, CBM-20A (article no. A59201)
- Gabi Star
- HERM flumo, HERM LB500
- Sedex85LT, Sedex90LT, Sedex100LT, Sedex LC (article no. A0754-x)
- CHIRALYSER-MP

Autosampler

AZURA® Column Thermostat CT 2.1 (article no. ATC00)

Fraction collectors

- Foxy® R1/R2 (article no. A59100/A59102/A591021)
- LABOCOL Vario-4000 (article no. A591022/ A591024)

Liquid Handler

- Liquid Handler LH 2.1 (article no. A5080)
- Liquid Handler LH 8.1

External pressure sensor (article no. AZG10-2)

Power plug Smartline series

Analytical Pumps 40P (article no. APC30xx)

Detectors

- UV Detector 40D
- UV Detector 50D

Smartline Degasser (article no. A5328)

Osmometers

This document is subject to technical changes.

Visit www.knauer.net for the latest version, ordering details and an overview of article numbers.